IFU SPECIFICATION

Specification # L8134-ZD
IFU Rev A
Document Rev 1
CO # 19-0668

IFU, TurboTemp EZ

Primary Supplier: Per current revision of Purchasing Master Table, as appropriate by product

Specification Requirements:

1. Color: Black and White
2. Languages: Arabic, Danish, Dutch, English, Finnish, French, German, Italian, Norwegian, Polish, Portuguese, Spanish, Swedish
3. Flat Size: 25.5" x 16.625"
4. Folded Size: 2.375" x 4.25"
5. Material: 27# white offset

Inspection:

1. Verify that the supplier is an approved supplier for this part number as stated above.
2. Inspect per IS-020.
3. Verify that the part number, lot number, and container quantity are affixed to each package, carton, container or bag.
4. Verify that the part number and lot number are detailed on the documentation included with the shipment.

Disposition:

Process in accordance with QAP-012.

CONFIDENTIAL: Zest Dental Solutions requires written permission for use or reproduction of this document in part or whole.

ENGLISH

TurboTemp EZ™ Self-Cure Provisional Composite

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

TurboTemp EZ™ is a self-cure, provisional crown-and-bridge composite that achieves long-lasting restorations with life-like esthetics. Developed as a multifunctional, highly-filled material, TurboTemp EZ is indicated for the fabrication of long- and short-term provisional restorations including bridges and long-span bridges, crowns, partial crowns, inlays/onlays, and veneers.

Featuring a high stress-at-yield advantage, accompanied by optimal flexibility, TurboTemp EZ is available in essential VITA shades that resemble the fluorescence of natural dentition, polishes beautifully, and can be trimmed, even at the margins, without breakage. With a predictable set time of less than 5 minutes, TurboTemp EZ offers reliable handling attributes that support the creation of visually pleasing, functional, and durable anterior and posterior restorations.

INDICATIONS FOR USE

- TurboTemp EZ is indicated for the fabrication of long- and short-term provisional restorations including veneers, inlays/onlays, crowns, partial crowns, bridges, and long-span bridges.

WARNINGS AND PRECAUTIONS

- Uncured TurboTemp EZ material can cause irritation, especially in people known to be allergic to methacrylates. Wash hands after handling the material.
- TurboTemp EZ is formulated to be used at room temperature. Shelf life at room temperature is stated on the cartridge.
- Turbo Temp EZ will bond to most dental adhesives and the oxygen-inhibited layer of resin-based restorative materials, making provisional removal for trimming difficult. Lightly lubricate with Vasoline or separating medium before fabrication.
- Do not cure TurboTemp EZ in proximity of eugenol-containing products. Do not let TurboTemp EZ come into contact with products containing eugenol. Eugenol can impair the hardening of TurboTemp EZ and cause discoloration.
- Do not expose TurboTemp EZ to elevated temperatures.
- Use a new mixing tip for every patient.

PRIOR TO APPLICATION

- Place the TurboTemp EZ cartridge into a 1:1 dispensing gun.
- Before each use, first bleed the cartridge until both the base and catalyst extrude out evenly.
- Then, attach a short blue mixing tip provided in the package and bleed the tip to ensure even mixing of the material.

WORK/REMOVE/SET TIMES

Work time (min:sec)	Removal from the tooth (min:sec)	Set time (min:sec)
0:45	2:30 – 3:00	5:00

Note: The times provided apply to a room temperature of 23° C and a normal relative air humidity of 50%. Different room temperature and humidity may result in variations from the above times. Higher temperatures may shorten and lower temperatures may increase these times.

APPLICATION

FABRICATION OF TEMPORARY RESTORATION

- Dispense the TurboTemp EZ material through the mixing tip into the impression. When extruding, keep the tip immersed in the material to avoid trapping air bubbles. Refer to the table above for information on working time.
 - Place the impression tray onto the prepared tooth and let it set for about 2:00 minutes from the start of mixing.
- Note:** The setting reaction can be monitored using an instrument on a pea-size material extruded on the mixing pad.
- Remove the temporary restoration from the impression tray during the elastic phase, within 2:30 – 3:00 minutes from the start of mixing.

Note: If oxygen-inhibited layer occurs, remove it with alcohol.

4. The material will fully set at about 5:00 minutes from the start of mixing.

FINISHING AND POLISHING OF TEMPORARY RESTORATION

- Remove the temporary restoration from the impression and complete the trimming using a dental bur or an instrument.
- Use the oxygen-inhibited layer using alcohol. The surface will appear matte, and will no longer be sticky when the restoration is fully cleaned.
- Achieve a final polish by using a flexible disk, rubber polisher, or goat's hair wheel until the restoration reaches a desired luster.

CEMENTATION OF TEMPORARY RESTORATION

- Attach the finished temporary restoration with a temporary luting cement.
- Note:** Eugenol cements can impair the hardening of TurboTemp EZ and cause discoloration.
- Check occlusion and adjust accordingly.

REPAIR AND CUSTOMIZATION OF TEMPORARY RESTORATIONS

Note: When contacting and repairing a freshly-made restoration, it does not need to be roughened, and a bonding agent is not required (proceed with step 3).

- Repair of an existing, or cleaned and polished restoration: slightly roughen the area to be modified using a grinder or sandblaster such as Zest's MicroEtch™.
- Apply a bonding agent such as Zest's Prelude One™ (to the prepared areas) according to the manufacturer's instructions for use.
- Apply additional TurboTemp EZ material, or a flowable material such as Zest's Accolade™ (to the prepared areas) according to the manufacturer's instructions for use.
- Finish and polish to complete the restoration.

STORAGE

Store at or below 77°F (25°C). Use at room temperature.

DEFINITIONS OF SYMBOLS

The following symbols may appear on the product packaging or labeling.

SYMBOL	TITLE	EXPLANATORY TEXT	STANDARD	REFERENCE
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC	ISO 15223-1	5.1.1
	European Community Authorized Representative	Indicates the authorized representative in the European Community	ISO 15223-1	5.1.2
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1	5.1.6
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1	5.1.5
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use	ISO 15223-1	5.4.3
	Use-by date	Indicates the date after which the medical device is not to be used	ISO 15223-1	5.1.4
	Upper Limit of Temperature	Indicates the temperature limits to which the medical device can be safely exposed	ISO 15223-1	5.3.6
	European Mark of Conformity	Indicates device is in conformance with Medical Device Directive 93/42/EEC	MDD 93/42/EEC	Annex XII
	Rx Only	Caution: Federal law restricts this device to sale by or on the order of a dentist	US Code of Federal Regulations, Title 21	801.15(c)(1) (i)(F)
	Quantity	Indicates the number of items within the package	N/A	N/A
	Medical device	Indicates the item is a medical device	DIS 15223-1 (2019)	5.7.7

	Danville Material 2875 136th Avenue East Carlsbad, CA 92010 USA www.zestdent.com	MDS5 GmbH Schifflingrain 41 30175 Hannover Germany	
---	--	--	---

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