<p> Turbulemp™ 3</p> <p> INSTRUCTIONS FOR USE</p> <p> PRODUCT DESCRIPTION</p> <p>Turbulemp™ 3 is a syringeable bi-acryl composite for chisade/provisional restorations. Turbulemp 3 is fast and accurate, especially when used in conjunction with a quality vinyl polysiloxane impression material such as Zest's First Quarter™, First Half™, or Start VPS™. Turbulemp 3 is available in 1/8, 1/4, 1/2, A3, A3.5, B1, and Bleach. All are delivered in 76g (2oz.) cartridges designed to fit on a 10:1 style automix gun. Ten waste save tips are included per refill kit.</p> <p> INDICATIONS FOR USE</p> <p>Turbulemp 3 is a temporary crown and bridge material designed as a syringeable acrylic-provisional composite for temporary crown and bridge placement.</p> <p> WARNINGS AND PRECAUTIONS</p> <ol style="list-style-type: none"> Turbulemp 3 contains methacrylate monomers which can cause allergic reactions in susceptible individuals. Avoid contact between uncured product and skin, oral soft tissues or eyes. Do not take internally. Use as directed. This product is intended for use by dental professionals only. Turbulemp 3 will adhesively bond to most dental adhesives and the air inhibited layer of resin-based restoratives, making provisional removal for trimming difficult. Do not store material in the proximity to eugenol-containing products, nor let the material come into contact with material containing eugenol. Eugenol can impair polymerization of the material and cause discoloration. Do not expose materials to elevated temperatures or intense light. Exposure to temperatures below 74°/723° will extend the setting time of Turbulemp 3. Set times are based on room temperature material. Refrigeration will delay set times. Avoid cross-contamination by using a new mixing tip for every patient. Do not use cartridge in-oral-ly. Do not use after expiration date. Turbulemp 3 MUST BE REMOVED FROM THE PATIENT PRIOR TO 30 DAYS. <p> TIMING</p> <p>0:00-0:40 – Insertion in the mouth 2:00-2:30 – Removal from the mouth 2:30-4:30 – Trimming/Finishing 5:00 – Final Hardness</p> <p> RECOMMENDED METHOD</p> <p> PRELIMINARY IMPRESSION: Prior to tooth preparation, place flexible vinyl polysiloxane (eg. Zest's First Quarter™ Monophase) on a TRIPLETRAY™ and make a closed bite impression. Slightly occlude material such as those the bite registration) must be avoided, as once removed, they will not be back well and will deteriorate. Always use an alternative, although less satisfactory.</p> <p> PREP AND FINAL IMPRESSION: Prepare subject tooth and complete final impression. To preclude bonding to Turbulemp 3, cover any composite buildup with a separating agent.</p> <p> IMMEDIATELY PRIOR TO USE: Remove cap and extrude a pea size quantity of material out of the cartridge. Eject slowly until a steady flow exudes from both compartments. Wipe off the end (without cross mixing) and install the mixing tip.</p> <p> FABRICATE TEMPORARY</p> <ol style="list-style-type: none"> CEMENTATION METHOD <ol style="list-style-type: none"> Discard the first pea size of mixed material and inject Turbulemp 3 into the prep areas of the preliminary impression (use care to avoid trapping air bubbles) and have patient close on the tray. Remove the tray when Turbulemp 3 has reached its elastic phase (2:00-2:30 minutes after injection). The provisional restoration will be retained in the tray and, for instance, Zest's Retraction Instrument. Move the instrument vertically around the margins while it is still in the tray using an amalgam carrier or #15 scalpel. Ensure composite is not locked into proximal undercuts. Immediately reet the tray and restoration in the mouth until completion of cure (approx. 5 minutes after injection). Remove the provisional restoration from the mouth and complete trimming and finishing with a diamond bur. Cement into place using a non-eugenol temporary cement, according to the manufacturer's instructions. Porosity can be filled with a flowable composite, according to the manufacturer's instructions. IMMEDIATELY AFTER USE: Remove cap and extrude a pea size quantity of material out of the cartridge. Eject slowly until a steady flow exudes from both compartments. Wipe off the end (without cross mixing) and install the mixing tip. SHRINK FIT METHOD <ol style="list-style-type: none"> Discard the first pea size of mixed material and inject Turbulemp 3 into the prep areas of the preliminary impression (use care to avoid trapping air bubbles) and have patient close on the tray. Remove the tray when Turbulemp 3 has reached its trimming/finishing phase (2:30-4:30 minutes after injection). Ideally the tray will come off the provisional restoration, leaving it firmly seated on the teeth. If not, immediately resect the restoration on the teeth. Trim off the flash using, for instance, Zest's Retraction Instrument. Move the instrument vertically to cleave the flash off the margins. Alternatively use a #15 scalpel blade to trim off the flash. Porosity can be filled with a flowable composite according to the manufacturer's instructions. For removal, the restoration will need to be split with a diamond bur and pieced off.
--

TURBTEMP 3 BRIDGE FABRICATION: Three units is the recommended maximum span. To add strength to the proximals of prosthetics, the conector area should be modified to fabricat, prior to taking the preliminary impression. In the posterior, both buccal and lingual can be modified. In the anterior, most of the modification should be done on the lingual to preserve esthetics. Use your preferred block-out material or soft wax.

ADDITIONAL REINFORCEMENT: Wet the fiber with Zest's E-Bond™ or Accolade™ flowable composite. Place the wetted fiber into the occlusal aspect of the preliminary impression. Using Turbulemp 3, infuse the fiber reinforcement and fill the remainder of the impression. Seat the filled preliminary impression in the mouth. Remove at approximately 2 minutes after injection. The reinforced provisional bridge will remain in the preliminary impression when it is removed from the mouth. Finish and cement as indicated in the recommended procedure.

HELPFUL HINTS




- When starting with a new cartridge: place cartridge in gun, remove cap, and extrude a small amount of material to ensure both sides are flowing. NOTE: Always bleed the cartridge before installing a new tip.
- Make sure to not move the mixing tip properly. The tip has different size bores and a notch to indicate proper orientation. Incorrectly mounting the tip can damage cartridge or cause wear remaining of material.
- A slightly gummy air inhibited layer can be seen on the hardened surface of the provisional restoration. This layer allows bubble and margin defects to be minimized by directly bonding with a flowable composite such as Zest's StartFlow™. The layer is easily removed with ethyl alcohol or polishing wheels/bushes.

STORAGE

Store at or below 77° F (25 °C). Do not freeze. Use at room temperature.

DEFINITIONS OF SYMBOLS

The following symbols may appear on the product packaging or labeling.

SYMBOL	TITLE	EXPLANATORY TEXT	STANDARD	REFERENCE
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC	ISO 15223-1	5.1.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community	ISO 15223-1	5.1.2
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1	5.1.6
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1	5.1.5
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use	ISO 15223-1	5.4.3
	Use-by date	Indicates the date after which the medical device is not to be used	ISO 15223-1	5.1.4
	Temperature Limitation	Indicates the temperature limits to which the medical device can be safely exposed	ISO 15223-1	5.3.6
	European Mark of Conformity	Indicates device is in conformance with Medical Directive 93/42/EEC	MD0 93/42/EEC	Annex XII
	Rx Only	Caution: Federal law restricts this device to sale by or on the order of a dentist	US Code of Federal Regulations, Title 21	801.15(c)(1) (OIF)
		MöSS GmbH 30175 Hannover/Germany		

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