

Flash Buster


Instructions for Use

Product Description

New composite burs for the removal of cement and stains; gentle on enamel, ceramics and soft tissue; self-sharpening, with a continuous-abrasive effect.

Available in different lengths, head shapes, and sizes.

To be used with a contra-angle (latch-type) handpiece, at a low speed (up to 10 000 r.p.m.), in the presence of water.



Indications for Use

- Removal of stains and colored coatings from the teeth surface, even difficult-to-reach areas, with no harm to the enamel.
- Cleaning of the enamel after ultrasonic scaling.
- Reduction of composite surplus.
- Removal of excess cement without denaturing the surrounding soft tissue.
- Trimming of temporary cement without altering preparations.
- In orthodontics: reduction of adhesive remnants from the tooth surface after removal of orthodontic brackets.
- In periodontics: after in-depth root scaling and cleaning to help support healing and flap reattachment.
- In implantology: easy removal of cement surplus from indirect restorations.

Precautions and Warnings

- Use these devices for indicated purposes only.
- Operate the devices up to the recommended speed only.





Cleaning and Sterilization

Warnings	Thoroughly dry the devices after all steps, if they will sit for more than 15 minutes prior to the next step. Follow all applicable solution or equipment manufacturer's instructions. All applicable solutions and equipment should meet national regulatory requirements. These devices are NOT supplied sterile. Clean and sterilize prior to use, as follows.
Processing limitations	None, other than wear, to be determined by the user.
Initial treatment at the point of use	Reprocess the devices as soon as possible after use. Do not allow material to dry. Keep the devices immersed in solution from the time

	they are used until they are cleaned – use enzymatic cleaner, cleaning/disinfecting, cleaning, or water (in order of preference).
Preparation before cleaning	Not applicable
Cleaning: Automated	1. Rinse the devices under cold, heavy running tap water. 2. Ultrasonically clean the devices in enzymatic cleaner for at least five minutes. Avoid the devices in contact with each other during cleaning (do not pile). 3. Rinse thoroughly with tap water.
Cleaning: Manual	Not applicable
Disinfection	Not applicable
Drying	Dry by wiping with absorbent, disposable wipes, dry heat, or filtered, pressurized air.
Maintenance, inspection and testing	Visually examine the devices for contamination. Magnification is recommended. Repeat the above steps for any devices that are still contaminated.
Packaging	Sterilization pouch
Sterilization	Steam sterilize, pre-vacuum, average pressure of 30psi (207 kPa), minimum of 4 minutes, 132°C, and a minimum dry time of 22 minutes.
Storage	Store the devices in the sterilized wrap.

Definition of Symbols

The following symbols may appear on the product packaging or labeling.

SYMBOL	TITLE	EXPLANATORY TEXT	STANDARD	REFERENCE
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC	EN 980 ISO 15223-1	5.12 5.1.1
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified	EN 980 ISO 15223-1	5.10 5.1.6
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified	EN 980 ISO 15223-1	5.4 5.1.5
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use	EN 980 ISO 15223-1	5.18 5.4.3
RxOnly	Rx Only	Caution: Federal law restricts this device to sale by or on the order of a dentist	US Code of Federal Regulations, Title 21	801.15(c)(1)(i)(F)



www.zestdent.com



Danville Materials
2875 Loker Avenue East
Carlsbad, CA 92010 USA
(1)760-743-7744