MicroCab™+ and MacroCab™+

INSTRUCTIONS FOR USE
PRODUCT DESCRIPTION
 The MicroCab™+ and MacroCab™+ dust cabinets are ideal for sandblast dust and handpiece grinding containment. It features a built in suction/filter system, as well as bright, LED lights for even illumination. The cabinets are compact and portable making it perfect for any lab or operator.

USE
 The cabinet is shipped with the power cord inside the unit. Prior to use, connect the power cord to the inlet on the back (see Figure 3) and to an electrical receptacle. The on/off switch under the integrated handle will operate the LED lights and suction/filter system.

When using the cabinet, only the fingertips with the MicroEtcher™ (Danville Materials) or handpiece need to be placed in the cabinet. Directing the abrasive spray towards the back of the cabinet will help contain the abrasive. It may tend to escape if the spray is directed towards a hand or operator.

PLACEMENT
 Keep the cabinet at least 6 inches (15 cm) away from any wall (back or sides) to ensure proper air flow and performance. The cabinet should be placed on a stable work surface or laboratory table.

Do not position the A/CDC adapter in such a way to block easy access to the power cord inlet (see Figure 3) or in any way that makes it difficult to disconnect the power cord.

Do not set the cabinet in a location where liquids or chemicals are used. Ensure that no liquids are spilled into the unit.

WINDOW INSTALLATION
 Four disposable windows are supplied with the original kit. Each window has a blue protective film covering it, which must be removed prior to use. The windows are installed by sliding the wide section up under the metal frame top (see Figure 1). By slightly bowing the window, the bottom edge is inserted into the front lip on the bottom of the window opening. Windows only need to be removed for replacement.


Figure 1

ELECTRICAL SHOCK HAZARD

For protection against electrical shock hazards, the equipment must be directly connected to an approved power supply with the power supply cord provided by the manufacturer. Do not use an ungrounded receptacle. Any break in the electrical ground path could be hazardous. If the power cord becomes cracked, frayed, broken, or otherwise damaged, please contact the manufacturer. Do not replace with a lower rated cord set.

CLEANING
 The cabinet is easily opened for cleaning. Cleaning is to be done regularly depending on frequency of use (see Figures 2 and 3).

- Unplug the power cord from the inlet on the back of the cabinet.
- Open the two latches on the back of the cabinet.
- Remove the window.
- Grasp the top handle, pull the filter and cabinet top straight up and out of the cabinet body.
- Dump all debris out of the cabinet. The filter can be removed and tapped against the sides of a trash pail. This allows powder trapped in the filter fibers to fall out. If the filter is removed from the filter assembly, reinstall it with the fins vertical.
- Reinstall the filter assembly by sliding it down into the cabinet. Ensure that the orange rubber of the filter slides behind the sealing fins that project into the cabinet from the side panels. There is also a channel under the front of the cabinet lid that captures the small rod running across the cabinet frame. This assures proper alignment.
- Reconnect the back cabinet latches and replace the window and power cord.


Figure 2


Figure 3

DISINFECTING
 Once the filter assembly is removed the interior of the cabinet can be sprayed with disinfectant and wiped down after abrasive/contaminants have been discarded.

Environmental Conditions
 The cabinet is designed for indoor operation in the following conditions:
 Altitude up to 2,000 m
 Temperature 5 °C to 40 °C
 Maximum relative humidity 80 % up to 31 °C decreasing linearly to 50 % at 40 °C
 Main supply voltage fluctuations of AC/DC adapter up to ±10% of nominal voltage
 Transient over voltages of AC/DC adapter up to the levels of OVERVOLTAGE CATEGORY II, POLLUTION DEGREE 2

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

TROUBLESHOOTING
 If you encounter any problems with operation of the cabinet, contact the manufacturer for help.


PRODUCT SPECIFICATIONS

MicroCab™+
 Size: 10.4" x 8.7" x 7.7" (depth x width x height)
 Weight: 5.7 lbs (2.6 kg)
 Power Supply AC Input: 100-240V, 1.5A, 50-60Hz
 Power Supply DC Output: 12.0V, 3.33A

MacroCab™+
 Size: 10.4" x 12.6" x 7.7" (depth x width x height)
 Weight: 8.4 lbs (3.8 kg)
 Power Supply AC Input: 100-240V, 1.5A, 50-60Hz
 Power Supply DC Output: 12.0V, 3.33A

DEFINITIONS OF SYMBOLS

SYMBOL	TITLE	EXPLANATORY TEXT	STANDARD	REFERENCE
	Manufacturer	Indicates the medical device manufacturer as defined in EU Directive 93/42/EEC	EN 980 ISO 15223-1	S.12 S.1.1
	Date of Manufacture	Indicates the date when the medical device was manufactured	EN 980 Annex A.7 Annex A.2	
	Catalogue Number	Indicates the manufacturer's catalogue number so that the medical device can be identified	EN 980 Annex A.7 Annex A.2	
	Batch Code	Indicates the manufacturer's batch code so that the batch or lot can be identified	EN 980 Annex A.7 Annex A.2	
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use or the instructions for use	EN 980 Annex A.7 Annex A.2	
	Caution, See Instructions for Use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself	EN 980 Annex A.7 Annex A.2	
	Carrying bag handles	To indicate the instructional safeguard to keep away from the carrying bag handles	IEC 60950-1	4.4.5
	Waste electrical and electronic equipment	Indicates that device requires separate collection at time of disposal should not be discarded with normal waste. The black bar indicates product placed on market after 2005.	WEEE Directive 2012/19/EU	Annex IX
	This way up	To indicate correct upright position of the transport packaging	ISO 780	13
	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully	ISO 15223-1	5.3.1
	European Mark of Conformity	Indicates device in conformance with Low Voltage Directive 2014/35/EU	LVD 2014/35/EU	n/a
	North American Safety Certification, "CS" indicates compliance for the "Canadian" standard "US" indicates conformance for the "US" standard.		OSHA/NRTL	n/a
	North American Safety Certification, "CS" indicates compliance for the "Canadian" standard "US" indicates conformance for the "US" standard.		OSHA/NRTL	n/a
	Danville Materials			

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