IFU SPECIFICATION

Specification # 89011-00
IFU Rev J
Document Rev 1
CO # 17-0712

Instructions for Use: Sure Etch

Product Name: INSTRUCTIONS FOR USE: Sure Etch

Primary Supplier: Robinson Printing

Specification Requirements:

1. Color: Black text on a white background.
2. English, Czech, German, Spanish, Finnish, French, Croatian, Hungarian, Italian, Dutch, Norwegian, Polish, Portuguese, Swedish languages
3. Flat size: 12.75" x 16.5"
4. Folded size: 4.25" x 2.75"

Inspection:

1. Verify that the supplier is an approved supplier for this part number as stated above.
2. Inspect per IS-020.
3. Verify that the Danville or Zest Dental Solutions part number, Danville or Zest Dental Solutions lot number, and container quantity are affixed to each package, carton, container or bag.
4. Verify that the Danville or Zest Dental Solutions part number and Danville or Zest Dental Solutions lot number are detailed on the documentation included with the shipment.

Disposition:

Process in accordance with QAP-012.

CONFIDENTIAL: Zest Dental Solutions requires written permission for use or reproduction of this document in part or whole.

ENGLISH

Sure Etch™ Instructions for Use

INDICATIONS FOR USE

Sure Etch™ is a 37% phosphoric acid etchant used to etch dentin or enamel prior to placement of a dental bonding agent.




ENAMEL AND DENTIN ETCHING

The appropriate etching time is determined by the manufacturer of the bonding agent being used. Therefore follow manufacturer’s instructions for etching. If etching instructions are not provided, proceed as follows:

- Clean the tooth surface with a pumice/water slurry or with aluminum oxide delivered by an air abrasion device such as PrepStart™ or PrepStart H₂O™ (Danville Materials).
- Rinse with water and dry with clean air.
- Apply Sure Etch to the tooth for 5-30 seconds, or as indicated by the adhesive manufacturer’s instructions for use:
 - Sure Etch Liquid is applied with an applicator brush, a straight bristled brush, or a sponge applicator.
 - Sure Etch Gel is applied directly out of the syringe after affixing a clean needle tip.
- Rinse well with water.
- Blow dry with oil-free air, or if appropriate for the bonding agent, blot dry or leave moist per the adhesive manufacturer’s instructions for use.
- Note:** When dry, etched enamel should have a frosted appearance. If the frosted appearance is absent, re-etch or check for etching barriers such as residual composite and take appropriate steps to remove it. Then re-etch.
- Proceed with the bonding agent application, per manufacturer’s instructions for use. It is important to avoid contaminating the etched tooth with saliva or other contaminants before the bonding agent is applied.
- After use, disinfect the Sure Etch container. If using Sure Etch Gel, discard the used tip and recap the syringe to prevent dehydration. If using Sure Etch Liquid, tightly recap the bottle to prevent evaporation.

DEFINITION OF SYMBOLS

The following symbols may appear on the product packaging or labeling.

SYMBOL	TITLE	EXPLANATORY TEXT	STANDARD	REFERENCE
	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC	EN 980 ISO 15223-1	5.12 5.1.1
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community	EN 980 ISO 15223-1	5.13 5.1.2
	Catalogue Number	Indicates the manufacturer’s catalogue number so that the medical device can be identified	EN 980 ISO 15223-1	5.10 5.1.6
	Batch Code	Indicates the manufacturer’s batch code so that the batch or lot can be identified	EN 980 ISO 15223-1	5.4 5.1.5
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use	EN 980 ISO 15223-1	5.18 5.4.3
	Use-by date	Indicates the date after which the medical device is not to be used	EN 980 ISO 15223-1	5.3 5.1.4
	Upper Limit of Temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed	EN 980 ISO 15223-1	5.17.2 5.3.6
	Keep Away from Sunlight/heat	Indicates a medical device that needs protection from light sources	EN 980 ISO 15223-1	5.20 5.3.2
	European Mark of Conformity	Indicates device is in conformance with Medical Device Directive 93/42/EEC	MDD 93/42/EEC	Annex XII
RxOnly	Rx Only	Caution: Federal law restricts this device to sale by or on the order of a dentist	US Code of Federal Regulations, Title 21	801.15(c)(1)(i)(F)

	Danville Materials 2875 Loker Avenue East Carlsbad, CA 92010 USA (1)760-743-7744 www.zestdent.com		MD55 GmbH Schiffgraben 41 30175 Hannover Germany		
		RxOnly		89011-00 REV J	
ZD ZEST DENTAL SOLUTIONS					

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