ENGLISH – Postupje i tehnike pripreme za Raymond Bertolotti, DDS, PhD. Za daljnje informacije kontaktirajte 5th Quarter Seminars na (510) 483-2411, FAKS (510) 652-8729.

Instructions for Use **RxOnly**
PRODUCT DESCRIPTION
First Quarter™ is a super hydrophilic, fast setting vinyl polysiloxane impression material intended to allow dentists to take accurate impressions of all types of tooth structures for various purposes, such as crown and bridge, implantology, etc. First Quarter is odorless, tasteless and immersible in disinfectants. It offers dimensional stability, tear resistance, and accuracy of impression. First Quarter materials are compatible with all other vinyl polysiloxane materials.

INDICATIONS
First Quarter is an additional-cure vinyl polysiloxane dental impression material that is used for all crown and bridges, edentulous and implant techniques.

First Quarter	Color	Room Temperature (23°C) min:sec	Setting Time Intraorally (35°C) min:sec
Light Body	Purple	0:30	1:45
Monophase (Chocomint Flavor)	Magenta	0:30	1:45
Monophase	Magenta	0:30	1:45
Heavy Body	Light Pink	0:30	1:45

Note: Work and set times may vary depending on storage conditions and temperature.

MIXING INSTRUCTIONS – CARTRIDGE

- Insert cartridge into an appropriate 1:1 mixing gun, remove twist-off cap, and extrude a small amount of material until even flow from both barrels is seen. Discard the dispensed material and wipe the cartridge end clean with tissue, avoiding cross contamination of the barrels.
- Attach the automix tip of the desired size. The material is now ready to use.
- After use, do not remove the automix tip. The used tip serves as a seal until removed for the next use.

CLOSED BITE IMPRESSIONS

The triple tray or closed bite impression is an efficient and accurate method to make an impression and establish bite registration.

- Tray Selection:
 - Anterior – Anterior Triple Tray™ (Premier)
 - Posterior – Sideless Triple Tray™ (Premier)
- Note:** Avoid rimmed posterior trays as they potentially induce distortions.

2. Technique:
An impression should be taken using two viscosities simultaneously; Monophase in the tray for dimensional stability and a vosc of Light Body for detail. Generally one person loads the tray while the second person syringes onto the tooth. The key to this procedure is to syringe Light Body only on clean, dry teeth, then blow with air until only a thin film remains. If a blank area remains, dry, syringe, and blow again, until only the thin film remains. Add Light Body to cover tooth, then seat tray.

Have patient close onto a tray of Monophase and guide patient into a CO closure. It is important to rehearse the proper closure beforehand.

Note: Putty should never be used for this procedure. It is too viscous, and induces elastic distortion.

It is critical that the Monophase be seated in the mouth before any elasticity develops. At additional working time is needed we recommend Start VPS™ in the normal set times. Heavy Body and Light Body Start VPS would be ideal.

SEPARATE FULL ARCH “PUTTY/WASH” IMPRESSIONS

(Use Light Body and Putty)
Creating accurate impressions using putty requires a dual set technique. Here, the putty is allowed to fully polymerize in the metal or plastic tray before the wash step.

- When using a custom tray made from a preliminary impression, use adhesive on the tray, following the adhesive manufacturer’s instructions for use. Light Body with the needle tip added to the mix tip is ideal.
- Before cutting the prep, make a putty impression, leaving room around the teeth for the wash. Leaving a space for the wash is achieved by simply placing a plastic film (such as a section of baggie or Reynolds Wrap) over the putty before seahing the tray. **Note:** Some plastic wraps will inhibit the set; test before use.
- Seat the tray with the putty, let polymerize, then remove tray and await prep.
- Use Light Body to take the final impression. Remove plastic film from the tray. Syringe Light Body onto clean dry teeth. Blow off with air until only a thin film remains. Repeat to cover any blank spots. The needle attachment for the small mixing tip is very handy for inlay, onlay and deep margins.
- Syringe Light Body into putty impression and seat.
- Remove after polymerization, wash and dry.

Note: Avoid simultaneous use; putty/wash set as putty is elastic and may cause distortion.

MONOPHASE IMPRESSION (USE MONOPHASE)

Single material impressions can be used where Light Body is not required for high flow. Monophase has a rapid set and fine texture, and is an ideal material to use for simple closed bite impressions as well as a preliminary for Turbo Temp™ (Zest Dental Solutions) temporary crown and bridge material.

- Syringe Monophase around clean, dry teeth. Syringe additional Monophase into sideless tray.
- Have the patient close until polymerized. Remove, wash and dry.

ADDITIONAL NOTES:

- First Quarter materials should be brought to room temperature prior to use. Exposure to prolonged temperatures above 77°F can be damaging. Store at room temperature.
- First Quarter gatives are compatible with all other vinyl polysiloxane materials.
- Some golves and particularly some powdered gloves can impart set. A sample test with a small amount is suggested. Keep putty jars closed when not in use.
- High viscosity materials used alone are not suitable for detailed impressions.
- Light Body impression materials used alone can fix excessively and may result in distortion.
- Procedures and techniques prepared courtesy of Raymond Bertolotti, DDS, PhD. For further information, please contact 5th Quarter Seminars at (510) 483-2411, FAX (510) 652-8729. www.adhesion.com

DEFINITIONS OF SYMBOLS

The following symbols may appear on the product packaging or labeling.

Symbol	Definition	Symbol	Definition
	Reference number		Lot number
	Use by		Manufacturer
RxOnly	Caution: Federal law restricts this device to sale by or on the order of a dentist.	EC REP	European Community Authorized Representative
	European Mark of Conformity		Consult instructions for use
	Upper limit of temperature		Keep away from sunlight

 <div>Damville Materials 2875 Lorch Avenue East Carlsbad, CA 92010 USA (1)760-743-7744</div>	EC REP	<div>MDSS GmbH Schiffgraben 41 30175 Hannover Germany</div>	
www.zestdent.com	RxOnly		89405 Rev G

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