

DANVILLE**IMPRESSION MATERIAL****Noginate™****INSTRUCTIONS**

Alginate Alternative: Noginate is an addition-curing elastomeric impression material intended to be used as an alternative to traditional alginate materials (e.g. for anatomical impressions.) Available in a 50ml cartridge.

CARTRIDGE PREPARATION

1. Insert cartridge into an appropriate 1:1 mixing gun, remove twist-off cap, and extrude a small amount of material until even flow from both barrels is seen. Discard the dispensed material and wipe the cartridge end clean with tissue, avoiding cross contamination of the barrels.
2. Attach an auto mix tip of desired size. The material is now ready for use.
3. After use, do not remove the automix tip. The used tip serves as a seal until removed for the next use.

INDICATIONS

Noginate is useful as alginate replacements for making impressions for:

- Preliminary impressions
- Anatomical models
- Fabricating temporary crowns and bridges
- Opposing dentition
- Fabricating simple removable prosthetic restorations

DANVILLE**IMPRESSION MATERIAL****Noginate™****INSTRUCTIONS**

Alginate Alternative: Noginate is an addition-curing elastomeric impression material intended to be used as an alternative to traditional alginate materials (e.g. for anatomical impressions.) Available in a 50ml cartridge.

CARTRIDGE PREPARATION

1. Insert cartridge into an appropriate 1:1 mixing gun, remove twist-off cap, and extrude a small amount of material until even flow from both barrels is seen. Discard the dispensed material and wipe the cartridge end clean with tissue, avoiding cross contamination of the barrels.
2. Attach an auto mix tip of desired size. The material is now ready for use.
3. After use, do not remove the automix tip. The used tip serves as a seal until removed for the next use.

INDICATIONS

Noginate is useful as alginate replacements for making impressions for:

- Preliminary impressions
- Anatomical models
- Fabricating temporary crowns and bridges
- Opposing dentition
- Fabricating simple removable prosthetic restorations

NOGINATE**INSTRUCTIONS**

- Orthodontic work
- Fabricating removable retainers and splints
- Case study models
- Producing models for the construction of trays for whitening solutions, bites and surgical guides.

ADDITIONAL NOTES:

- Noginate materials should be brought to room temperature prior to use. Exposure to prolonged temperatures above 77°F can be damaging. Store at room temperature.
- Noginate material is compatible with all other vinyl polysiloxane materials.

 Danville Materials
a Zest Anchors, LLC company
2875 Loker Avenue East
Carlsbad, CA 92010 USA
(1)760-743-7744

www.zestdent.com

ZD ZEST DENTAL
SOLUTIONS

91635 REV D

NOGINATE**INSTRUCTIONS**

- Orthodontic work
- Fabricating removable retainers and splints
- Case study models
- Producing models for the construction of trays for whitening solutions, bites and surgical guides.

ADDITIONAL NOTES:

- Noginate materials should be brought to room temperature prior to use. Exposure to prolonged temperatures above 77°F can be damaging. Store at room temperature.
- Noginate material is compatible with all other vinyl polysiloxane materials.

 Danville Materials
a Zest Anchors, LLC company
2875 Loker Avenue East
Carlsbad, CA 92010 USA
(1)760-743-7744

www.zestdent.com

ZD ZEST DENTAL
SOLUTIONS

91635 REV D