

# VPS Impression Materials

## Ideal Materials For Every Technique & Speed

Highest quality and value priced. Developed with the help of leading clinicians to make consistent, reliable, results easy! Wash materials are super-hydrophilic and give a high-gloss surface and high-tear strength. Precise setting time and snap-set. Placement in mouth triggers an accelerated set. Formulated from the highest quality polymers and fillers.

Available in six color coded viscosities for easy differentiation: Super Flow, Light Body, Monophase, Heavy Flex, Heavy Stiff Bite, Clear Bite & Putty.

The material is free of taste and is suitable for a variety of impression-taking techniques.

- Putty and Super Flow Light Body formulations can be employed for the standard putty-wash technique; the smooth putty resists sticking to fingers and gloves.
- Heavy Stiff Bite is the perfect bite registration material and provides a high degree of accuracy.
- Heavy Flex can be used for closed bite impressions (using Triple Tray impression trays) in combination with Super Flow Light Body for closed-bite impressions.
- The very fast-setting Monophase material also is formulated for Triple Tray impressions and is particularly suitable for those that do not involve hard-to-reach areas (which generally require a super flow viscosity).



## 3 SET TIMES



### First Quarter Fast Set Impression Material

A super fast impression material great for smaller, Triple Tray type impressions.



### First Half Medium Set Impression Material

Ideal for those that want speed, but need more time than is offered by the First Quarter series. New Orthodontic Monophase is ideal for single viscosity full arch impressions. Available in 50ml cartridges or 380ml cartridges for Pentamix.










### Star VPS Regular Set Impression Material

offers a variety of viscosities in the standard set times, leaving plenty of working time for multiple unit impressions!

First Quarter	Working Time  1/2 Minute s	Total Time To Set  3/4 Minute s	Durometer	Color	2 Pack4	Pack	Econo 24 Pack	Super Econo 96 Pack
Light Body			48		89277	89272	89462	90731
Monophase (Chocomint Flavor)			62		90901	90900	90902	90903
Monophase			62		89386	89387	80008-01	90729
Heavy Body			65		90913	90809	90810	90811

First Half	Working Time  1 1/4 Minute s	Total Time To Set  2 1/4 Minute s	Durometer	Color	2 Pack	4 Pack	Econo 24 Pack	Super Econo 96 Pack
Light Body			48		90467	90468	90469	90733
Light Body			48		91241	91242	91243	91244
Monophase Ortho			62		91122	91123	91124	91125
Heavy Flex			63		90648	90649	90650	90732
Putty-FS			73		93133	(2 400gm Tubs Base & Catalyst)		

Star VPS	Working Time	Total Time To Set	Durometer	Color	2 Pack	4 Pack	Econo 24 Pack	Super Econo 96 Pack		
Ultra Light			48		N/A	80060-01	80061-01	N/A		
Light Body			48		80021-01	80011-01	80009-01	90727		
Heavy Flex			63		80023-01	80013-01	80007-01	90728		
Ultra Heavy			65N		/A	80062-01	80065-01	N/A		
Heavy Stiff Bite			20 Seconds	50 Seconds	90		80024-01	80014-01	80006-01	90730
Clear Bite			25 Seconds 2	Minutes	63		80025-01	80015-01	80037-01	N/A
Putty			90 Seconds	4 1/2 Minutes	65		80016	(2 400gm Tubs Base & Catalyst)		
AlphaBite			30 Seconds	60 Seconds	82		93352	N/A	93354N	/A
Noginate	90 Seconds	2 1/2 Minutes	45		93553	N/A	93554N	/A		