

TurboTemp™ 3

Temporary Crown & Bridge Material



ADVANTAGES

- 91% clinical rating
- Easy to trim
- Very good esthetics
- Stain resistant
- Polishes easily
- Easy to extrude

TurboTemp3 is a further improvement with better physical properties than even Luxatemp®* & Protemp®*, yet is still about 30% less costly.

Turbotemp 3 is in a 10:1 cartridge which aids in the superior physical properties.

TurboTemp3 produces provisional restorations that have very good reproduction of detail and polish to a smooth finish. Danville will send a free slider to convert the current users with a 4:1 gun to a 10:1 gun (this is a value of \$19.99)



TT3: Starter Kit: 76 gm Cartridge, 10 Mix Tips & Dispensing Gun					
A-1	A-2	A-3	A-3.5	B-1	Bleach
93319	93320	93379	93321	93322	93592

TT3: Refills: 76 gm Cartridge & 10 Mix Tips					
A-1	A-2	A-3	A-3.5	B-1	Bleach
93323	93324	93369	93325	93326	93591

*Not a Danville Material trademark