

NEW


CHAIRSIDE®
Impression Materials

INTRODUCING
CHAIRSIDE® VPS IMPRESSION MATERIALS



CHAIRSIDE® VPS IMPRESSION MATERIALS

IDEAL MATERIALS FOR YOUR IMPLANT CASES

Here we grow again! ZEST Anchors is expanding the CHAIRSIDE Product Portfolio with a new line of VPS Impression Materials ideal for your implant cases. Introducing **CHAIRSIDE® VPS Impression Materials** - High quality, value conscious materials developed with leading clinicians for ease of use while producing consistent, reliable results.

- Formulated from the **highest quality** polymers and fillers
- Exceptional **accuracy** and dimensional **stability** with high tear resistance
- Clinician friendly working times and patient friendly **quick** set times

CHAIRSIDE Impression Material

- Wash material is super-hydrophilic and gives a high-gloss surface and high-tear strength
- Ideal for full arch or triple tray impressions
- Available in 3 color-coded viscosities for easy differentiation while accommodating for a variety of impression techniques



CHAIRSIDE Bite Registration

- Highly accurate, stable and non-brittle
- Quick intra-oral setting time
- Easily trimmable

Here is what your peers are saying about
CHAIRSIDE VPS IMPRESSION MATERIALS

I found the CHAIRSIDE VPS Impression Material products from ZEST Anchors to produce consistent, reliable results, in addition the cost is very competitive resulting in an overall value for my practice."

Charles Mastrovich, DDS, Escondido, California, USA

IDEAL MATERIALS FOR A VARIETY OF IMPRESSION TECHNIQUES

Material Properties for CHAIRSIDE® Impression Materials

	TEAR STRENGTH	CONSISTENCY	ELASTIC RECOVERY	SHORE A HARDENESS (15 MINUTES)	TENSILE STRENGTH / ELONGATION
Light	3.8 MPa	40mm	99%	48	3.8 MPa / 224.9%
Monophase	4.5 MPa	36mm	99%	62	3.8 MPa / 164.6%
Heavy	4.2 MPa	34mm	99%	63	4.3 MPa / 158.3%
Bite Registration	6.7 MPa	30mm	99%	90	6.9 MPa / 26.7%

Ordering Information

ITEM NUMBER	DESCRIPTION	MATERIAL COLOR	WORKING TIME	SETTING TIME	PRICE
09569	CHAIRSIDE Impression Material, Light Body (includes (2) 50ml cartridges, 6 mixing tips, 6 curved tips)		1 Minute	2 ½ Minutes	\$40.00
09570	CHAIRSIDE Impression Material, Medium Body (Monophase) (includes (2) 50ml cartridges, 6 mixing tips)		1 Minute	2 ½ Minutes	\$40.00
09571	CHAIRSIDE Impression Material, Heavy Body (includes (2) 50ml cartridges, 6 mixing tips)		1 Minute	2 ½ Minutes	\$40.00
09572	CHAIRSIDE Bite Registration Material (includes (2) 50ml cartridges, 6 mixing tips)		20 Seconds	50 Seconds	\$40.00





CHAIRSIDE[®]
Impression Materials



Try **CHAIRSIDE VPS Impression Materials** today
by calling 800.262.2310 or visiting www.zestanchors.com



ZEST ANCHORS
2061 WINERIDGE PLACE
ESCONDIDO, CA 92029
USA

TEL: 760.743.7744 OR 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTANCHORS.COM
WWW.ZESTANCHORS.COM

DISTRIBUTED BY: