

NEW


CHAIRSIDE®
Block Out Material

**IDEAL FOR
BLOCKING OUT UNDERCUTS
WHILE PROCESSING
OVERDENTURES CHAIR SIDE**




ZEST
ANCHORS



BLOCK OR AVOID UNDERCUTS USING CHAIRSIDE® BLOCK OUT MATERIAL

ZEST Anchors continually listens to clinicians about how the company can design, manufacture and source products that make implant-retained overdenture treatment simpler and more efficient, as well as creating a comfortable and satisfying procedure for patients. The latest advancement in the CHAIRSIDE Product Portfolio supports that mission - **Introducing CHAIRSIDE Block Out Material** from ZEST Anchors, a **light-cured protective material that removes undercuts while processing overdentures chair side.**




PATIENT COMFORT

CHAIRSIDE Block Out Material is an acrylic resin that remains cool and comfortable during light curing, preventing discomfort and harm to the tissue.

EASY HANDLING

Viscous material stays put using syringe delivery. Contrasting blue color is easy to see, flexible and simple to remove.



MULTIPLE USES

CHAIRSIDE Block Out is a superior material for blocking out undercuts for LOCATOR® cases - it is also ideal to:

- Protect tissue while bleaching, sandblasting or using air abrasion
- Seal around rubber dams
- Repair torn rubber dams
- Maintain reservoir space for bleaching trays

ITEM NUMBER	DESCRIPTION	PRICE
09574	<ul style="list-style-type: none"> • Four (4) 1.2mL Syringes • Twenty (20) Mixing Tips 	\$34.99



ZEST ANCHORS
2061 WINERIDGE PLACE
ESCONDIDO, CA 92029
USA

TEL: 760.743.7744 OR 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTANCHORS.COM
WWW.ZESTANCHORS.COM