

  
**LOCATOR<sup>®</sup>**  
OVERDENTURE IMPLANT SYSTEM

# **LOCATOR<sup>®</sup>** **OVERDENTURE** **IMPLANT SYSTEM**

PRODUCT CATALOG



**ZD ZEST DENTAL**  
SOLUTIONS<sup>®</sup>

# THE LOCATOR® OVERDENTURE IMPLANT SYSTEM

**\$179<sup>99</sup>**

ALL-INCLUSIVE PACKAGE CONTAINS

- 1 Implant
- 1 LOCATOR Abutment
- 1 Processing Pack

**\$79<sup>99</sup>**


IMPLANT ONLY PACKAGE




SAVE OVER \$70 WITH ALL-IN-ONE!

# THE LOCATOR® OVERDENTURE IMPLANT SYSTEM

## PROCESSING PACK



Denture Housing



Blue Standard Range Insert  
Low Retention





Pink Standard Range Insert  
Medium Retention



Red Extended Range Insert  
Low Retention



Block-Out Spacer



2.4mm			2.9mm			3.5mm		
Part	Length	Cuff	Part	Length	Cuff	Part	Length	Cuff
07450	10mm	2.5mm	07460	10mm	2.5mm	07501-02	8mm	2.5mm
07451	12mm	2.5mm	07461	12mm	2.5mm	07502-02	10mm	2.5mm
07452	14mm	2.5mm	07462	14mm	2.5mm	07503-02	12mm	2.5mm
07440	10mm	3mm	07443	10mm	3mm	07504-02	14mm	2.5mm
07441	12mm	3mm	07444	12mm	3mm	07501-03	8mm	3mm
07442	14mm	3mm	07445	14mm	3mm	07502-03	10mm	3mm
07455	10mm	4mm	07465	10mm	4mm	07503-03	12mm	3mm
07456	12mm	4mm	07466	12mm	4mm	07504-03	14mm	3mm
07457	14mm	4mm	07467	14mm	4mm	07501-04	8mm	4mm
07432	10mm	5mm	07435	10mm	5mm	07502-04	10mm	4mm
07433	12mm	5mm	07436	12mm	5mm	07503-04	12mm	4mm
07434	14mm	5mm	07437	14mm	5mm	07504-04	14mm	4mm
07381	10mm	Implant Only	07386	10mm	Implant Only	07501-05	8mm	5mm
07382	12mm	Implant Only	07387	12mm	Implant Only	07502-05	10mm	5mm
07383	14mm	Implant Only	07388	14mm	Implant Only	07503-05	12mm	5mm
						07504-05	14mm	5mm
						07501	8mm	Implant Only
						07502	10mm	Implant Only
						07503	12mm	Implant Only
						07504	14mm	Implant Only

3.9mm			4.4mm			4.9mm		
Part	Length	Cuff	Part	Length	Cuff	Part	Length	Cuff
07505-02	8mm	2.5mm	07509-02	8mm	2.5mm	07513-02	8mm	2.5mm
07506-02	10mm	2.5mm	07510-02	10mm	2.5mm	07514-02	10mm	2.5mm
07507-02	12mm	2.5mm	07511-02	12mm	2.5mm	07515-02	12mm	2.5mm
07508-02	14mm	2.5mm	07512-02	14mm	2.5mm	07516-02	14mm	2.5mm
07505-03	8mm	3mm	07509-03	8mm	3mm	07513-03	8mm	3mm
07506-03	10mm	3mm	07510-03	10mm	3mm	07514-03	10mm	3mm
07507-03	12mm	3mm	07511-03	12mm	3mm	07515-03	12mm	3mm
07508-03	14mm	3mm	07512-03	14mm	3mm	07516-03	14mm	3mm
07505-04	8mm	4mm	07509-04	8mm	4mm	07513-04	8mm	4mm
07506-04	10mm	4mm	07510-04	10mm	4mm	07514-04	10mm	4mm
07507-04	12mm	4mm	07511-04	12mm	4mm	07515-04	12mm	4mm
07508-04	14mm	4mm	07512-04	14mm	4mm	07516-04	14mm	4mm
07505-05	8mm	5mm	07509-05	8mm	5mm	07513-05	8mm	5mm
07506-05	10mm	5mm	07510-05	10mm	5mm	07514-05	10mm	5mm
07507-05	12mm	5mm	07511-05	12mm	5mm	07515-05	12mm	5mm
07508-05	14mm	5mm	07512-05	14mm	5mm	07516-05	14mm	5mm
07505	8mm	Implant Only	07509	8mm	Implant Only	07513	8mm	Implant Only
07506	10mm	Implant Only	07510	10mm	Implant Only	07514	10mm	Implant Only
07507	12mm	Implant Only	07511	12mm	Implant Only	07515	12mm	Implant Only
07508	14mm	Implant Only	07512	14mm	Implant Only	07516	14mm	Implant Only


# SURGICAL KIT AND INSTRUMENTATION



Part	Product	Description	Price
07357		LOCATOR® Implant Latch Driver	\$69.99
07360		LOCATOR® Implant Driver, Short	\$69.99
07361		LOCATOR® Implant Driver, Long	\$69.99
07362		Torque Indicating Wrench	\$403.99
07463		Torque Wrench Insert	\$99.99
08926		LOCATOR® Square Driver Torque Wrench Insert	\$55.99
08393		LOCATOR® Core Tool	\$94.99
08390		LOCATOR® Core Tool Abutment Driver (Gold End)	\$44.99
08397		LOCATOR® Core Tool Abutment Driver (Tip End)	\$38.99
08394		LOCATOR® Abutment Holder Sleeve (4 Pack)	\$24.99
07374		Drill Extender	\$29.99
07365		Direction Indicator	\$18.99
07373		LOCATOR® Implant Rotary Tissue Punch, 3.7mmD	\$104.99
07534		LOCATOR® Implant Rotary Tissue Punch, 5.8mmD	\$104.99
07537		Starter Drill, 1.5mm	\$59.99
07538		● Drill, 1.6mm	\$59.99
07539		● Drill, 2.1mm	\$59.99
07540		● Drill, 2.4mm	\$59.99
07541		● Drill, 2.6mm	\$59.99
07542		● Drill, 2.9mm	\$59.99
07543		○ Drill, 3.2mm	\$59.99
07544		● Drill, 3.6mm	\$59.99
07545		● Drill, 4.1mm	\$59.99
07546		● Drill, 4.6mm	\$59.99



# SURGICAL KIT AND INSTRUMENTATION

Part	Product	Description	Price
07533		Narrow Ridge Drill Stops 4, 6, 8, 10, 12, 14mm	\$69.99
07536		Standard Ridge Drill Stops 4, 6, 8, 10, 12, 14mm	\$69.99
07556		Cortical Drill 3.3mm	\$59.99
07557		Cortical Drill 3.7mm	\$59.99
07558		Cortical Drill 4.3mm	\$59.99
07559		Cortical Drill 4.8mm	\$59.99
07560	Standard Surgical Tray		\$529.99
07561	LOCATOR® Implant Surgical Kit		\$1,299.99
07564	LOCATOR® Implant Premium Surgical Kit		\$1,699.99
09120		LOCATOR® Implant Patient Education Model	\$189.99









# LOCATOR® ACCESSORY COMPONENTS





## MEASUREMENT COMPONENTS

Part	Product	Description	Price
09530		Angle Measurement Guide	\$26.99
08517		Parallel Post (4 Pack)	\$46.99

## PROCESSING COMPONENTS

Part	Product	Description	Price
08505		Impression Coping (4 Pack)	\$78.99
08505-20		Impression Coping (20 Pack)	\$349.99
08530		Abutment Analog (4mm) (4 Pack)	\$78.99
08530-20		Abutment Analog (4mm) (20 Pack)	\$354.99
07339		Healing Cap, 3mm	\$34.99
07340		Healing Cap, 4mm	\$34.99
08515		Black Processing Replacement Male (4 Pack)	\$29.99
08515-20		Black Processing Replacement Male (20 Pack)	\$135.99
08514		Block-Out Spacer (20 Pack)	\$20.99
08569		Processing Spacer (4 Pack)	\$32.99
08569-20		Processing Spacer (20 Pack)	\$157.99


## TORQUE WRENCH & DRIVERS

Part	Product	Description	Price
09020		30Ncm Torque Wrench, Thumb, Knob and 15mm Locator Insert Driver	\$264.99
04391	20Ncm Torque Wrench, Thumb, Knob and 15mm Locator Insert Driver		\$264.99
08913		Torque Wrench Insert Driver - 23mm total length (Short) (Latch Type)	\$55.99
08914		Torque Wrench Insert Driver - 29mm total length (Long) (Latch Type)	\$55.99
08926		LOCATOR® Square Driver Torque Wrench Insert Driver (15mm)	\$55.99
08927		LOCATOR® Square Driver Torque Wrench Insert Driver (21mm)	\$55.99
04936		.050 Hex Insert Driver for 30Ncm Torque Wrench (Square Drive Connection) (21mm)	\$55.99

# LOCATOR® ACCESSORY COMPONENTS


## STANDARD INSERTS

Dual retention to maximize stability and pivoting action that accommodates up to 20° of total divergence between two implants.











## EXTENDED RANGE INSERTS

Pivoting action accommodates up to 40° of total divergence between two implants.



## LOCATOR RETENTION INSERTS & DENTURE ATTACHMENT HOUSINGS

Part	Product	Description	Price
08510		Denture Attachment Housing Assembly (4 Pack)	\$69.99
08510-10		Denture Attachment Housing Assembly (10 Pack)	\$159.99
08529		Blue: Low Retention Standard Range Inserts (4 Pack)	\$29.99
08529-20		Blue: Low Retention Standard Range Inserts (20 Pack)	\$135.99
08527		Pink: Medium Retention Standard Range Inserts (4 Pack)	\$29.99
08527-20		Pink: Medium Retention Standard Range Inserts (20 Pack)	\$135.99
08524		Clear: High Retention Standard Range Inserts (4 Pack)	\$29.99
08524-20		Clear: High Retention Standard Range Inserts (20 Pack)	\$135.99
08548		Red: Low Retention Extended Range Inserts (4 Pack)	\$29.99
08548-20		Red: Low Retention Extended Range Inserts (20 Pack)	\$135.99
08915		Orange: Medium Retention Extended Range Inserts (4 Pack)	\$29.99
08915-20		Orange: Medium Retention Extended Range Inserts (20 Pack)	\$135.99
08547		Green: High Retention Extended Range Inserts (4 Pack)	\$29.99
08547-20		Green: High Retention Extended Range Inserts (20 Pack)	\$135.99
08558		Gray: Zero Retention Extended Range Inserts (4 Pack)	\$29.99
08558-20		Gray: Zero Retention Extended Range Inserts (20 Pack)	\$135.99

## ATTACHMENTS FOR LOCATOR IMPLANTS

Part	Description	Price
07351	2.5mm cuff height	\$125.99
07318	3mm cuff height	\$125.99
07352	4mm cuff height	\$125.99
07319	5mm cuff height	\$125.99
07320	6mm cuff height	\$125.99

# THE ZEST DENTAL SOLUTIONS® 10 YEAR WARRANTY



## THE ZEST DENTAL SOLUTIONS NARROW AND STANDARD RIDGE IMPLANT WARRANTY: LOCATOR® OVERDENTURE IMPLANT SYSTEM ("LODI").

Zest Dental Solutions ("Zest") is committed to providing quality products and dedicated to gathering feedback about its products. Zest actively collects and reviews the feedback of users of our products in compliance with regulatory reporting requirements and to better help us understand market expectations and validate our products' performance. The collection method Zest utilizes for such feedback is the Zest Product Experience Report ("PER") form. The PER form is to be completed from information provided by the attending clinician to share their Zest product experience.

Pursuant to the Zest implant warranty, Zest will replace covered LODI implants such qualifying covered implant is returned (as such terms are defined below). Upon a request for replacement under the Zest warranty, Zest will send a replacement product once Zest receives the returned product and completed PER form and confirms that the returned product is covered under the Zest warranty.

ONLY DIRECT ZEST CUSTOMERS MAY MAKE A WARRANTY CLAIM. In order to make a claim under the Zest warranty please return the items mentioned below in protective packaging and send these via a shipping method which enables the package to be tracked:

Printed copy of PER form completed from information provided by the attending clinician.

Explanted product(s) in sterile condition (NON-STERILE products will not qualify for replacement)

If a product failure or loss of integration has occurred, please additionally send relevant radiographs (these will not be returned unless specifically requested, please send copies)

Send shipment to: Zest Dental Solutions  
ATTN: Customer Service (US customers) or Wholesale Distribution  
(Distributors/OEM Partners)  
2875 Loker Ave East Carlsbad, CA 92010

If using the electronic form, please include the Tracking Number for the returned product package.

Download Zest NDI Product Experience Report (PDF) AT  
WWW.ZESTDENT.COM

Zest Dental Solutions ("Zest") IMPLANT WARRANTY (Valid as of May 1, 2015)

**1. Warranty beneficiary and scope:** Zest Dental Solutions ("Zest") hereby warrants to the direct Zest customer purchasing the from Zest ("Customer") that the Zest LOCATOR Overdenture Implant ("LODI"), when implanted according to the respective Zest LODI Technique Manual and other written instructions provided by Zest by a clinician (the "User"), will be free from any loss or lack of integration, fracture or other structural failure for the period of 10 years ("Warranty Period") from the time of treatment by the User (collectively, the "Zest Warranty"). This warranty only applies to the Customer. Third parties, particularly patients, are not covered by the Zest warranty and have no rights hereunder. Customers' sole remedy and Zest's sole liability under this Zest warranty is the replacement of the LODI implant by Zest as set forth herein. The Zest implant warranty only covers the replacement of the LODI implant and not any associated costs or expenses, including, but not limited to, chair time, laboratory fees and any other associated treatment.

**2. Zest Warranty conditions:** In the event that any request for warranty service for an implant is made by a Customer under the Zest warranty during the Warranty Period, Zest will replace such implant with the same or substantially equivalent product subject to the terms and conditions herein. The replacement product will be sent upon receipt of the returned product and the completed PER form and confirmation that the product is covered under the Zest warranty. To qualify for coverage under the Zest warranty, the claim must be made within the applicable Warranty Period by a Customer and all conditions below must be met. Once these requirements are satisfied, Zest will send the replacement part(s). The following REQUIRED conditions must be met and documented in order for coverage of a returned implant under

this Zest implant warranty:

**1.** The LODI was used exclusively with all components, connections, attachments and other technology provided by Zest and not in combination with any other manufacturer's products or technology;

**2.** The LODI is returned in sterile condition (or disinfected if delivered as such);

**3.** The LODI was implanted by a User and inspected and maintained in full compliance with the respective Zest LODI Technique Manual valid at the time of treatment and all other Zest written instructions as well as recognized dental procedures, during and after the treatment;

**4.** The patient had good oral hygiene which was monitored and documented by the User

**5.** The LODI was not subjected to damage caused by misuse, misapplication, accident, trauma or any other damage caused by external factors or the User, the patient or a third party;

**6.** A completed and signed PER was completed and submitted no later than 10 days after the complaint is made. Customer is responsible to ensure that the PER is completed and submitted from information provided by the applicable User. Details of the incident are imperative to determine whether vigilance reporting is required to regulatory authorities.

**3. Limits and limitations:** This Zest implant warranty is the only guarantee provided by Zest and shall apply in addition to any warranty rights conferred under any written sales agreement executed by Zest. The User remains free to claim rights against his supplier. EXCEPT AS SET FORTH HEREIN AND IN A WRITTEN SALES AGREEMENT EXECUTED BY ZEST, ZEST HEREBY DISCLAIMS ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, WITH RESPECT TO LODI/SNDI OR ANY OTHER ZEST PRODUCTS, SERVICES OR INFORMATION, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY OF FITNESS FOR PURPOSE OR MERCHANTABILITY, OR NON-INFRINGEMENT OF THIRD PARTY RIGHTS, OR ANY WARRANTY BY COURSE OF PERFORMANCE OR OTHERWISE.

ZEST SHALL NOT BE LIABLE TO THE CUSTOMER, THE USER, THE PATIENT, OR ANY THIRD PARTY, FOR LOST EARNINGS OR PROFITS, DIRECT OR INDIRECT DAMAGES AS WELL AS COLLATERAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES, DIRECTLY OR INDIRECTLY RELATED TO LODI OR ANY OTHER ZEST PRODUCTS, SERVICES OR INFORMATION.

**4. Zest Warranty territory:** This Zest implant warranty applies worldwide to LODI implants sold by Zest, a Zest affiliated company or an official distributor of Zest on or after the validity date stated above.

**5. Modification or termination:** Zest may modify or terminate this Zest implant warranty at any time in whole or in part. Changes to, or the termination of the Zest implant warranty, will not affect the warranty given for LODI installed prior to the date of the change or termination.

### LOCATOR OVERDENTURE IMPLANT (LODI) PRODUCT EXCHANGE POLICY

Zest Dental Solutions understands that customers may need to adjust inventories of LOCATOR Overdenture Implant System (LODI) Products in order to achieve the correct mix of sizes to treat their patients. Zest Dental Solutions will waive normal restocking fees (1:1 or greater) for exchanges during the first six months following the original purchase. Packaging cannot be written on or in any way adulterated. Shipping will still be the responsibility of the customer requesting the exchange.

# RETURN POLICY AND WARRANTY

## ORDERING

Orders are accepted by internet, email, mail, fax, or phone. Regular business hours are 5:00am to 5:00pm PST (Pacific Standard Time). Domestic orders received by 1:00pm PST can be shipped by ground the same day. Orders received by 3:00pm PST can be shipped by express. Institutional orders require a purchase order number.

To place a new order, please choose one of the following ways:

1. Online: [www.zestdent.com](http://www.zestdent.com)
2. E-mail: [sales@zestdent.com](mailto:sales@zestdent.com)
3. Phone: 1-800-262-2310
4. Fax: 1-800-487-1357

## PRICING

Zest Dental Solutions makes every effort to maintain a competitive pricing structure. Pricing is subject to routine review and change without prior notice. All prices listed are in U.S. dollars.

## TERMS AND BILLING

- A. Credit Card:** Payment at the time of order with VISA, AMERICAN EXPRESS, DISCOVER or MASTERCARD is necessary for customers not approved for open account billing or special financing offers.
- B. C.O.D.:** Delivery can also be made by C.O.D. freight collect.
- C. Open Account:** An open account can be established by completing a credit application and receiving approval by Zest Dental Solutions. The terms of orders purchased under open accounts are net 30 days.
- D. Export Orders:** Individual orders shipped outside the U.S. require payment in advance or a letter of credit.
- E. Past Due Accounts:** Past due balances will be subject to a 1.5% finance charge per month, amount equal to 18% per annum.

## SHIPPING CHARGES

Merchandise is shipped prepaid by Zest Dental Solutions with the cost added to the invoice, or freight collect with C.O.D. fee included in the case of C.O.D. shipments. Delivery options include Fed Ex ground, 3-day, 2-day, or overnight rush service.

## LIMITED WARRANTY

Zest Dental Solutions provides a limited warranty for its products, to the original purchaser, to be free from defects in workmanship and materials under normal use and service, for a period of one year from the date of purchase with the exception of LOCATOR and LOCATOR R-Tx abutments which carry a three year limited warranty (when purchased directly from Zest Dental Solutions) and a two year limited warranty when purchased from one of our partners. Zest Dental Solutions will, at its option, substitute the returned product that proves defective within the warranty period, with a similar product, free of charge.

Zest Dental Solutions continually strives to improve its products, and therefore, reserves the right to improve, modify, or discontinue products at any time without notice or incurring obligations. Purchaser assumes all risks and liability resulting from the use of Zest Dental Solutions Products, whether used separately or in combination with other products not of Zest Dental Solutions manufacture.

## RETURN POLICY

Please Observe the Following Guidelines:

- A.** Authorization for returns must be received from Zest Dental Solutions prior to the return of any product. A Return Material Authorization (RMA) number will be provided for all returns.
- B.** Shipping charges must be prepaid by the customer to accept a return shipment.
- C.** Returned products are subject to a \$20 restocking fee for orders under \$400. For orders \$400 to \$1,000, a \$40 restocking fee will be applied. For orders exceeding \$1,000, please contact your Customer Service Representative for the restocking fee.
- D.** Returned merchandise will be accepted within 90 days of purchase if product is in saleable condition (in its original unopened package and not marred by any added writing or over-labeling).
- E.** Returns will not be accepted after 90 days of purchase.
- F.** Non-returnable items include:
  1. Merchandise retained beyond expiration date noted on the package.
  2. Packages with broken seals or missing parts
  3. Used, damaged, or obsolete products will not be accepted for return.
- G.** LOCATOR Implant returns will be accepted in a 1:1 exchange for a period of 6 months with no restocking fee.



VISIT OUR WEBSITE AT  
**WWW.ZESTDENT.COM**  
TO PLACE ORDERS ONLINE  
24 HOURS A DAY 7 DAYS A WEEK



**ZEST DENTAL SOLUTIONS**  
2875 LOKER AVENUE EAST  
CARLSBAD, CA 92010  
USA

TEL: 1.855.866.LODI (5934)  
FAX: 760.743.7975  
EMAIL: SALES@ZESTDENT.COM  
**WWW.ZESTDENT.COM**

DISTRIBUTED BY: