



CHAIRSIDE®

Product Portfolio

INTRODUCING

CHAIRSIDE® IMPRESSION AND BITE REGISTRATION MATERIALS



CHAIRSIDE® VPS IMPRESSION AND BITE REGISTRATION MATERIALS IDEAL MATERIALS FOR ALL IMPLANT CASES

CHAIRSIDE® VPS Impression and Bite Registration Materials are additional-cure vinyl polysiloxane (VPS) materials developed to capture accurate impressions for all implant cases, in a predictable and consistent manner.

Available in medium set-times and color-coded viscosities with light, medium, and heavy-body formulations, these materials allow practitioners to efficiently address a variety of clinical situations including challenging, long-span cases and edentulous treatments.

CHAIRSIDE Impression Materials

- Wash material is super-hydrophilic and gives a high-gloss surface and high-tear strength
- Ideal for full-arch or triple-tray impressions
- Available in 3 color-coded viscosities for easy differentiation while accommodating a variety of impression techniques

CHAIRSIDE Bite Registration

- Highly accurate, stable and non-brittle
- Quick intra-oral setting time
- Easily trimmable



Features & Benefits:

- Exceptional accuracy
- High tear strength and resistance
- Excellent dimensional stability
- Hydrophilic and immersible in disinfectants

Here is what your peers are saying about
CHAIRSIDE VPS MATERIALS

I found the CHAIRSIDE VPS Impression Material products from Zest Dental Solutions to produce consistent, reliable results, in addition the cost is very competitive resulting in an overall value for my practice.”





Charles Mastrovich, DDS, Escondido, California, USA


IDEAL MATERIALS FOR A VARIETY OF IMPRESSION TECHNIQUES

Material Properties of CHAIRSIDE® Materials

CHAIRSIDE Impression and Bite Registration Materials	COLOR	WORKING TIME At Room Temperature (73°F/23°C) min:sec	SETTING TIME Intraorally (95°F/35°C) min:sec	DUROMETER READING
Light Body	Magenta	2:15 minutes	4 minutes	43
Medium Body	Dark Blue			48
Heavy Body	Yellow			63
Bite Registration	Light Blue	20 seconds	60 seconds	85

Product Ordering Information

ITEM NUMBER	DESCRIPTION	MATERIAL COLOR	PRICE
09614	CHAIRSIDE Impression Material, Light Body (includes (2) 50ml cartridges, 6 mixing tips, 6 curved tips)		\$40.00
09615	CHAIRSIDE Impression Material, Medium Body (includes (2) 50ml cartridges, 6 mixing tips)		\$40.00
09616	CHAIRSIDE Impression Material, Heavy Body (includes (2) 50ml cartridges, 6 mixing tips)		\$40.00
09572	CHAIRSIDE Bite Registration Material (includes (2) 50ml cartridges, 6 mixing tips)		\$40.00





CHAIRSIDE®

Product Portfolio



Try CHAIRSIDE® VPS Impression and Bite Registration Materials today by calling 800.827.7940 or visiting www.zestdent.com



ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 760.743.7744 OR 800.827.7940
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM

DISTRIBUTED BY: