

LOCATOR R-Tx[®] REMOVABLE ATTACHMENT SYSTEM



U.S. RETAIL PRODUCT CATALOG AND PRICE LIST*

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

System Overview and Abutment Selection Criteria

Each LOCATOR R-Tx Attachment System is delivered to you in a convenient, all-in-one double-ended vial that separately holds the Abutment and Processing Components.


Order for your favorite implant systems using this product ordering and reference guide.

LOCATOR R-Tx
Abutment



\$155 PRICE
INCLUDES:


- 1 - Abutment
- 1 - Denture Attachment Housing with Processing Insert
- 1 - Zero (**Gray**) Retention Insert
- 1 - Low (**Blue**) Retention Insert
- 1 - Medium (**Pink**) Retention Insert
- 1 - High (**Clear**) Retention Insert
- 1 - Block-Out Spacer



Denture Attachment Housing
with Black Processing Insert



Zero (**Gray**)
Retention Insert



Low (**Blue**)
Retention Insert



Medium (**Pink**)
Retention Insert




High (**Clear**)
Retention Insert



Block-out
Spacer

SELECTING THE PROPER CUFF HEIGHT FOR THE LOCATOR R-Tx ABUTMENT



1. Using a periodontal probe, measure from the top of the implant to the highest point of the gingiva.
2. Select the cuff height of the LOCATOR R-Tx Abutment that corresponds to that measurement. For example, if the measurement is 2mm, select the 2mm cuff. If the measurement is 2.5mm, select the 3mm cuff. An additional 2mm of Abutment height will extend above the gingiva to accommodate for the Denture Attachment Housing and Retention Insert.

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information

LOCATOR R-Tx ABUTMENTS ARE IMPLANT SPECIFIC

TO ORDER, YOU WILL NEED
THE FOLLOWING INFORMATION:

1. Name of the implant system
2. Implant diameter
3. Connection type and platform diameter
4. Tissue measurement (at the highest point)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

ACE SURGICAL

BONE LEVEL CONNECTION (Compatible with infinity Octagon Implants)

Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06

TISSUE LEVEL CONNECTION (Compatible with infinity Octagon Implants)

Regular Platform (RP)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
-----------------------	--	----------	----------	----------	----------	----------	----------

INTERNAL HEX CONNECTION (Compatible with infinity Internal HEX Implants)

3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06

TRI-CHANNEL CONNECTION (Compatible with infinity TRI-CAM Implants)

3.5mm Diameter Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Diameter Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Diameter Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

EXTERNAL HEX CONNECTION (Compatible with infinity EXTERNAL HEX Implants)

3.5mm Diameter Platform		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Diameter Platform		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm Diameter Platform (Platform Switch Only)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06

BIOHORIZONS

INTERNAL HEX CONNECTION (Compatible with Tapered Plus, Tapered Internal Implant, Tapered Tissue Level, Single Stage and Internal Parallel Wall Implants)

3.0mm Diameter Platform		30600-01	30600-02	30600-03	30600-04	30600-05	30600-06
3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

BIOMET 3i NOTE: BIOMET 3i Abutments feature a .048 hex drive

CERTAIN® CONNECTION (Compatible with 3i T3®, NanoTite™ and OSSEOTITE® Implants)

3.4mm Diameter Certain Platform		30300-01	30300-02	30300-03	30300-04	30300-05	30300-06
4.1mm Diameter Certain Platform		30301-01	30301-02	30301-03	30301-04	30301-05	30301-06

EXTERNAL HEX CONNECTION (Compatible with 3i T3®, NanoTite™ and OSSEOTITE® Implants)

3.4mm Diameter External Hex Platform		30305-01	30305-02	30305-03	30305-04	30305-05	30305-06
4.1mm Diameter External Hex Platform		30306-01 (1.2mm)	30306-02	30306-03	30306-04	30306-05	30306-06

BLUE SKY BIO

BONE LEVEL CONNECTION (Compatible with BIO Quattro Implants)

Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06

TISSUE LEVEL CONNECTION (Compatible with BIO One Stage Implants)

Regular Platform (RP)		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06
Wide Platform (WP)		30907-01	30907-02	30907-03	30907-04	30907-05	30907-06

CAMLOG (Call BioHorizons for availability - 888.246.8338)

DENTSPLY

ANKYLOS CONICAL CONNECTION (Compatible with Ankylos® C/X Implants)

Ankylos Implant			30715-02	30715-03	30715-04	30715-05	30715-06
-----------------	--	--	----------	----------	----------	----------	----------

ASTRA TECH EV CONNECTION (Compatible with Astra Tech EV Implants)

3.6mm Diameter Implant		30700-01	30700-02	30700-03	30700-04	30700-05	
4.2mm Diameter Implant		30701-01	30701-02	30701-03	30701-04	30701-05	
4.8mm Diameter Implant		30702-01	30702-02	30702-03	30702-04	30702-05	

ASTRA TECH TX CONNECTION (Compatible with Astra Tech TX Implants)

Astra TX (Aqua)		30706-01	30706-02	30706-03	30706-04	30706-05	
Astra TX (Lilac)		30707-01	30707-02	30707-03	30707-04	30707-05	

XIVE CONNECTION (Compatible with XiVE® Implants)

3.4mm Diameter Implant		30711-01	30711-02	30711-03	30711-04	30711-05	
3.8mm Diameter Implant		30712-01	30712-02	30712-03	30712-04	30712-05	
4.5mm Diameter Implant		30713-01	30713-02	30713-03	30713-04	30713-05	
5.5mm Diameter Implant		30714-01	30714-02	30714-03	30714-04	30714-05	

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

NEW

HIOSSEN

EXTERNAL HEX CONNECTION (Compatible with USII SA, USIII SA, Implants)

Mini Platform (Yellow)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
Regular Platform (RP) (Green)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Wide Platform (WP) (Blue)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06

INTERNAL HEX CONNECTION (Compatible with ET Implants)

Mini Platform (Yellow)		31600-01	31600-02	31600-03	31600-04	31600-05	31600-06
Standard Platform (SP) (Green)		31601-01	31601-02	31601-03	31601-04	31601-05	31601-06

NEW

IMPLANT DIRECT

INTERNAL HEX CONNECTION (Compatible with Legacy Implants)

3.5mm Diameter Platform (Green)	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform (Purple)	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform (Yellow)		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06

TRI-CHANNEL CONNECTION (Compatible with RePlant, RePlus and ReActive Implants)

3.5mm Diameter Platform (Purple)	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Diameter Platform (Yellow)	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Diameter Platform (Blue)		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

CONICAL CONNECTION (Compatible with Interactive Implants)

3.0mm Diameter Platform (Purple)		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
3.4mm Diameter Platform (Yellow)		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06

INTERNAL HEX CONNECTION (Compatible with Legacy2, Legacy3 and Legacy4 Implants)

3.0mm Diameter Platform		31500-01 (0.75mm)	31500-02	31500-03	31500-04	31500-05	31500-06
-------------------------	--	----------------------	----------	----------	----------	----------	----------

INTERNAL HEX CONNECTION (Compatible with Swish Implants)

4.8mm Diameter Platform		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06
-------------------------	--	----------	----------	----------	----------	----------	----------

NEW

KEYSTONE DENTAL

TILOBE CONNECTION (Compatible with Genesis, PrimaConnex, and MAX - TL 5.0/5.5/6.5mm Platform Implants)

3.5/3.8mm Diameter (Green) Platform		31200-01	31200-02	31200-03	31200-04	31200-05	31200-06
4.1/4.5mm Diameter (Blue) Platform		31201-01	31201-02	31201-03	31201-04	31201-05	31201-06
5.0/5.5/6.5mm Diameter (Purple) Platform		31202-01	31202-02	31202-03	31202-04	31202-05	31202-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

NOBEL BIOCARE

CONICAL CONNECTION (Compatible with NobelActive®, NobelParallel™ CC and NobelReplace® CC Implants)

Narrow Platform (NP)		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
Regular Platform (RP)		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06
Wide Platform (WP)		30508-01	30508-02	30508-03	30508-04	30508-05	30508-06

EXTERNAL CONNECTION (Compatible with Brånemark and NobelSpeedy® Groovy Implants)

Narrow Platform (NP)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
Regular Platform (RP)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
Wide Platform (WP)		30512-01	30512-02	30512-03	30512-04	30512-05	

TRI-CHANNEL CONNECTION (Compatible with NobelReplace®, ReplaceSelect™ and NobelSpeedy® Replace Implants)

Narrow Platform (NP)	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
Regular Platform (RP)	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
Wide Platform (WP)		30502-01	30502-02	30502-03	30502-04	30502-05	

OCO BIOMEDICAL

INTERNAL HEX CONNECTION (Compatible with Engage® Implants)

NEW

3.0mm Diameter Platform		31500-01 (0.75mm)	31500-02	31500-03	31500-04	31500-05	31500-06
3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06

SOUTHERN IMPLANTS

TRI-NEX CONNECTION (Compatible with Tri-Nex™, Tri-Nex Co-Axis® and Tri-MAX 5.0mm Platform Implants)

3.5mm Tri-Nex Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Tri-Nex Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Tri-Nex Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

NEW

STRAUMANN

BONE LEVEL (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

NC Platform		30401-01	30401-02	30401-03	30401-04	30401-05	30401-06
RC Platform		30402-01	30402-02	30402-03	30402-04	30402-05	30402-06

TISSUE LEVEL (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

NNC Platform		30405-01	30405-02	30405-03	30405-04	30405-05	30405-06
RN Platform		30406-01	30406-02	30406-03	30406-04	30406-05	30406-06
WN Platform		30407-01 (1.5mm)	30407-02	30407-03	30407-04	30407-05	30407-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

ZIMMER

INTERNAL HEX CONNECTION	(Compatible with Tapered Screw-Vent®, Screw-Vent®, Trabecular Metal™ Implants and AdVent® 5.7mm Platform Implant)						
3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06

LOCATOR R-Tx[®] REMOVABLE ATTACHMENT SYSTEM

Accessory Components

PART #	DESCRIPTION		PRICE
PROCESSING COMPONENTS			
30001-01	Zero Retention Insert (Gray), (QTY. 4)		\$25.00
30002-01	Low Retention Insert (Blue), (QTY. 4)		\$25.00
30003-01	Medium Retention Insert (Pink), (QTY. 4)		\$25.00
30004-01	High Retention Insert (Clear), (QTY. 4)		\$25.00
30012-01	Processing Insert (QTY. 4)		\$25.00
08514	Block-Out Spacer (QTY. 20)		\$20.00
30013-01	Denture Attachment Processing Assembly (QTY. 4)		\$60.00
30014-01	3.35mm Abutment Analog (QTY. 4)		\$75.00
30015-01	4mm Abutment Analog (QTY. 4)		\$75.00
30016-01	5mm Abutment Analog (QTY. 4)		\$75.00
30017-01	Impression Coping (QTY. 4)		\$75.00
30018-01	Processing Spacer (QTY. 4)		\$34.00
TOOLS			
30021-01	LOCATOR R-Tx Retention Insert Tool (Dual Sided)		\$55.00
08007	.050"/1.25mm Hex Driver with thumb knob (Shaft 6.4mm)		\$86.00
08008	.050"/1.25mm Hex, Latchlock Driver (Shaft 12.7mm)		\$53.00
04936	.050"/1.25mm Hex, Torque Wrench Insert (Shaft 6.4mm)		\$53.00

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM


Accessory Components (Continued)

PART #	DESCRIPTION	PRICE
PATIENT EDUCATION		
09665	LOCATOR R-Tx Patient Education Model	\$180.00

09665




LOCATOR R-Tx[®]
REMOVABLE ATTACHMENT SYSTEM



VISIT **ZESTDENT.COM**
24 HOURS A DAY | 7 DAYS A WEEK

TO PLACE ORDERS ONLINE OR TO FIND
THE LATEST SUPPORT TOOLS FOR
LOCATOR R-Tx[®]

ZD ZEST DENTAL SOLUTIONS[™]
ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM

DISTRIBUTED BY: