



LOCATOR F-Tx[®] **FIXED** ATTACHMENT SYSTEM

U.S. RETAIL PRODUCT CATALOG
AND PRICE LIST*



*Products and pricing listed in this guide may not be available in all countries.


LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

System Overview and Abutment Selection Criteria

Each LOCATOR F-Tx Abutment is delivered in convenient all-in-one packaging that includes all the sterile components needed for each implant.

Order for your favorite implant systems using this product ordering and reference guide.

LOCATOR F-Tx
Abutment




\$249 PRICE
INCLUDES:

- 1 - Abutment
- 1 - Denture Attachment Housing with Pre-inserted Processing Ball
- 1 - Low (**Blue**) Retention Ball
- 1 - Medium (**Tan**) Retention Ball
- 1 - High (**Green**) Retention Ball
- 1 - Additional Processing Ball
- 2 - Block-Out Spacers (Packaged separately)




SELECTING THE PROPER CUFF HEIGHT FOR THE LOCATOR F-Tx ABUTMENT



1. Using a periodontal probe, measure from the top of the implant to the highest point of the gingiva.
2. The cuff size should match the tissue height, as this will place the prosthetic finish line right at the gingival height and maximizes the retention surface of the Denture Attachment Housing in the prosthesis. The cuff margin is at the widest portion of the sphere (midpoint of the Abutment sphere). For subgingival finish line, choose a shorter cuff size.

LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

Abutment Ordering Information



LOCATOR F-Tx ABUTMENTS ARE IMPLANT SPECIFIC

TO ORDER, YOU WILL NEED THE FOLLOWING INFORMATION:

- 1.** Name of the implant system
- 2.** Implant diameter
- 3.** Connection type and platform diameter
- 4.** Cuff height selection criteria

IMPLANT COMPANY	LOCATOR F-Tx ATTACHMENT SYSTEM (Includes Abutment)					
DESCRIPTION	CUFF HEIGHT					
	1.5mm	2mm	3mm	4mm	5mm	6mm

ACE SURGICAL

TRI-CHANNEL (Compatible with infinity TRI-CAM Implants)						
3.5mm Diameter Platform	10500-01	10500-02	10500-03	10500-04	10500-05	10500-06
4.3mm Diameter Platform	10501-01	10501-02	10501-03	10501-04	10501-05	10501-06
5.0mm Diameter Platform	10502-01	10502-02	10502-03	10502-04	10502-05	10502-06
EXTERNAL HEX CONNECTION (Compatible with infinity EXTERNAL HEX Implants)						
3.5mm Diameter Platform		10510-02	10510-03	10510-04	10510-05	10510-06
4.1mm Diameter Platform		10511-02	10511-03	10511-04	10511-05	10511-06
5.0mm Diameter Platform (Platform Switch Only)		10511-02	10511-03	10511-04	10511-05	10511-06
INTERNAL HEX CONNECTION (Compatible with infinity Internal HEX Implants)						
3.5mm Diameter Platform	10200-01	10200-02	10200-03	10200-04	10200-05	10200-06
4.5mm Diameter Platform	10201-01	10201-02	10201-03	10201-04	10201-05	10201-06
BONE LEVEL CONNECTION (Compatible with infinity Octagon Implants)						
Narrow Platform (NP)	10900-01	10900-02	10900-03	10900-04	10900-05	10900-06
Regular Platform (RP)	10901-01	10901-02	10901-03	10901-04	10901-05	10901-06
TISSUE LEVEL CONNECTION (Compatible with infinity Octagon Implants)						
Regular Platform (RP)		11302-02	11302-03	11302-04	11302-05	11302-06

BIOHORIZONS

INTERNAL HEX CONNECTION (Compatible with Tapered Plus, Tapered Internal Implant, Tapered Tissue Level, Single Stage and Internal Parallel Wall Implants)						
3.5mm Internal Hex Platform	10200-01	10200-02	10200-03	10200-04	10200-05	10200-06
4.5mm Internal Hex Platform	10201-01	10201-02	10201-03	10201-04	10201-05	10201-06

LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR F-Tx ATTACHMENT SYSTEM (Includes Abutment)					
DESCRIPTION	CUFF HEIGHT					
	1.5mm	2mm	3mm	4mm	5mm	6mm

BIOMET 3i

CERTAIN® INTERNAL CONNECTION (Compatible with 3i T3®, NanoTite® and OSSEOTITE® Implants)						
3.4mm Diameter CERTAIN Platform	10300-01	10300-02	10300-03	10300-04	10300-05	10300-06
4.1mm Diameter CERTAIN Platform	10301-01	10301-02	10301-03	10301-04	10301-05	10301-06
EXTERNAL HEX CONNECTION (Compatible with 3i T3®, NanoTite® and OSSEOTITE® Implants)						
4.1mm Diameter External Hex Platform		10511-02	10511-03	10511-04	10511-05	10511-06

BLUE SKY BIO

BONE LEVEL CONNECTION (Compatible with BIO Quattro Implants)						
Narrow Platform (NP)	10900-01	10900-02	10900-03	10900-04	10900-05	10900-06
Regular Platform (RP)	10901-01	10901-02	10901-03	10901-04	10901-05	10901-06
TISSUE LEVEL CONNECTION (Compatible with One Stage Implants)						
Regular Platform (RP)		10906-02	10906-03	10906-04	10906-05	10906-06

NEW CAMLOG

CAMLOG CONNECTION (Compatible with CAMLOG Implants)						
3.3mm Diameter Platform	10800-01	10800-02	10800-03	10800-04	10800-05	10800-06
3.8mm Diameter Platform	10801-01	10801-02	10801-03	10801-04	10801-05	10801-06
4.3mm Diameter Platform	10802-01	10802-02	10802-03	10802-04	10802-05	10802-06
5.0mm Diameter Platform	10803-01	10803-02	10803-03	10803-04	10803-05	10803-06
iSy CONNECTION (Compatible with iSy Implant)						
3.5mm Diameter Platform	10810-01	10810-02	10810-03	10810-04	10810-05	10810-06

DENTSPLY

ANKLYOS® CONNECTION (Compatible with Ankylos® C/X Implants)						
Ankylos Implant Platform		10715-02	10715-03	10715-04	10715-05	10715-06
ASTRA TECH EV CONNECTION (Compatible with Astra Tech EV Implants)						
3.6mm Conical EV	10700-01	10700-02	10700-03	10700-04	10700-05	10700-06
4.2mm Conical EV	10701-01	10701-02	10701-03	10701-04	10701-05	10701-06

LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR F-Tx ATTACHMENT SYSTEM (Includes Abutment)					
DESCRIPTION	CUFF HEIGHT					
	1.5mm	2mm	3mm	4mm	5mm	6mm

DENTSPLY

ASTRA TECH TX CONNECTION (Compatible with Astra Tech TX Implants)						
ASTRA TX (AQUA)	10706-01	10706-02	10706-03	10706-04	10706-05	10706-06
ASTRA TX (LILAC)	10707-01	10707-02	10707-03	10707-04	10707-05	10707-06

KEYSTONE DENTAL

TiLOBE CONNECTION (Compatible with Genesis, PrimaConnex and MAX TL 5.0/5.5/6.5mm Platform Implants)						
3.5/3.8mm Diameter(Green) Platform	11200-01	11200-02	11200-03	11200-04	11200-05	11200-06
4.1/4.5mm Diameter(Blue) Platform	11201-01	11201-02	11201-03	11201-04	11201-05	11201-06
5.0/5.5/6.5mm Diameter(Purple) Platform	11202-01	11202-02	11202-03	11202-04	11202-05	11202-06

NEODENT

MORSE TAPER CONNECTION (Compatible with CM Implants)						
CM Implant Platform		10715-02	10715-03	10715-04	10715-05	10715-06

NOBEL BIOCARE

CONICAL CONNECTION (Compatible with NobelActive®, NobelParallel™ CC and NobelReplace® CC Implants)						
Narrow Platform (NP)	10506-01	10506-02	10506-03	10506-04	10506-05	10506-06
Regular Platform (RP)	10507-01	10507-02	10507-03	10507-04	10507-05	10507-06
Wide Platform (WP)	10508-01	10508-02	10508-03	10508-04	10508-05	10508-06

EXTERNAL HEX CONNECTION (Compatible with Brånemark and NobelSpeedy™ Groovy Implants)						
Narrow Platform (NP)		10510-02	10510-03	10510-04	10510-05	10510-06
Regular Platform (RP)		10511-02	10511-03	10511-04	10511-05	10511-06
Wide Platform (WP)		10512-02	10512-03	10512-04	10512-05	10512-06

TRI-CHANNEL CONNECTION (Compatible with NobelReplace®, Replace Select™ and NobelSpeedy™ Replace Implants)						
Narrow Platform (NP)	10500-01	10500-02	10500-03	10500-04	10500-05	10500-06
Regular Platform (RP)	10501-01	10501-02	10501-03	10501-04	10501-05	10501-06
Wide Platform (WP)	10502-01	10502-02	10502-03	10502-04	10502-05	10502-06

LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR F-Tx ATTACHMENT SYSTEM (Includes Abutment)					
DESCRIPTION	CUFF HEIGHT					
	1.5mm	2mm	3mm	4mm	5mm	6mm

SOUTHERN IMPLANTS

TRI-NEX CONNECTION (Compatible with Tri-Nex™, Tri-Nex Co-Axis® and Tri-MAX 5.0mm and 6.0mm Platform Implants)

3.5mm Tri-Nex Platform	10500-01	10500-02	10500-03	10500-04	10500-05	10500-06
4.3mm Tri-Nex Platform	10501-01	10501-02	10501-03	10501-04	10501-05	10501-06
5.0mm Tri-Nex Platform	10502-01	10502-02	10502-03	10502-04	10502-05	10502-06
6.0mm Tri-Nex Platform	11713-01	11713-02	11713-03	11713-04	11713-05	11713-06

EXTERNAL HEX CONNECTION (Compatible with External Hex, External Hex Co-Axis and External Hex MAX Implants)

4.0mm Diameter Platform		11703-02	11703-03	11703-04	11703-05	11703-06
5.0mm Diameter Platform		11704-02	11704-03	11704-04	11704-05	11704-06
5.7mm Diameter Platform		11705-02	11705-03	11705-04	11705-05	11705-06

STRAUMANN

BONE LEVEL CONNECTION (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

NC Platform	10401-01	10401-02	10401-03	10401-04	10401-05	10401-06
RC Platform	10402-01	10402-02	10402-03	10402-04	10402-05	10402-06

TISSUE LEVEL CONNECTION (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

NNC Platform		10405-02	10405-03	10405-04	10405-05	10405-06
RN Platform		10406-02	10406-03	10406-04	10406-05	10406-06
WN Platform		10407-02 (2.5mm)	10407-03	10407-04	10407-05	10407-06












ZIMMER BIOMET

INTERNAL HEX CONNECTION (Compatible with Tapered Screw-Vent®, Screw-Vent®, Trabecular Metal™ and AdVent® 5.7mm Platform Implants)

3.5mm Diameter Platform	10200-01	10200-02	10200-03	10200-04	10200-05	10200-06
4.5mm Diameter Platform	10201-01	10201-02	10201-03	10201-04	10201-05	10201-06
5.7mm Diameter Platform	10202-01	10202-02	10202-03	10202-04	10202-05	10202-06

LOCATOR F-Tx[®] FIXED ATTACHMENT SYSTEM

Ancillary Components

PART #	DESCRIPTION	PRICE	
PROCESSING COMPONENTS			
10041-01	Denture Attachment Housing (Qty. 2)		\$80.00
10032-01	Black Processing Ball (Qty. 2)		\$35.00
10033-01	Low Retention Ball (Blue), (Qty. 2)		\$35.00
10034-01	Medium Retention Ball (Tan), (Qty. 2)		\$35.00
10035-01	High Retention Ball (Green), (Qty. 2)		\$35.00
10036-01	Healing Cap (Qty. 2)		\$40.00
10037-01	Abutment Analog (Qty. 2)		\$40.00
10040-01	Waxing Cap (Qty. 2)		\$40.00
10043-01	Impression Coping (Qty. 2)		\$80.00
10030-01	Processing Cap (Qty. 2)		\$40.00
10039-01	Polishing Cap (Qty. 2)		\$40.00

LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

Instrumentation and Tools (Continued)

PART #	DESCRIPTION	PRICE
DRIVERS		
10028-01	LOCATOR F-Tx Abutment Driver, Manual (used with square drive connections)	\$55.00
10029-01	LOCATOR F-Tx Abutment Driver, Latch (used with latch-lock connections)	\$55.00
10042-01	LOCATOR F-Tx Retention Ball Hex Driver	\$55.00

10100










LOCATOR F-Tx SEATING AND REMOVAL TOOL

10100	LOCATOR F-Tx Seating and Removal Tool (includes Tool, Seating Tip, Removal Tip with Threaded Adapter, Wire/Lever Wrench and Tip Wrench).	\$239.00
10094	LOCATOR F-Tx Tip Wrench	\$12.99
10095	LOCATOR F-Tx Seating Tip	\$19.99

LOCATOR F-Tx® FIXED ATTACHMENT SYSTEM

Instrumentation and Tools (Continued)

PART #	DESCRIPTION	PRICE
10096	LOCATOR F-Tx Removal Tip with Threaded Adapter 	\$99.99
10097	LOCATOR F-Tx Wire Replacement 	\$19.99
10098	LOCATOR F-Tx Threaded Adapter 	\$19.99
10099	LOCATOR F-Tx Wire and Lever Wrench 	\$19.99
PROSTHESIS REMOVAL TOOL BAR AND LOOP		
10023	LOCATOR F-Tx Removal Loop 	\$30.00
10024	LOCATOR F-Tx Removal Tool Bar 	\$99.00
PATIENT EDUCATION		
09673	LOCATOR F-Tx Patient Education Model 	\$360.00



VISIT
ZESTDENT.COM
24 HOURS A DAY | 7 DAYS A WEEK

TO PLACE ORDERS ONLINE OR TO FIND THE LATEST
SUPPORT TOOLS FOR LOCATOR F-Tx®

ZD ZEST DENTAL SOLUTIONS
ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM

DISTRIBUTED BY: