


NEW



CHAIRSIDE[®]
Denture Removal Tool

**GIVE YOUR
PATIENTS THE GIFT
OF EASY DENTURE
REMOVAL**



**ZD ZEST DENTAL
SOLUTIONS™**

ZEST | DANVILLE MATERIALS | PERIOSCOPY



CHAIRSIDE®

Denture Removal Tool


DO YOU HAVE PATIENTS OR PRIMARY CAREGIVERS THAT STRUGGLE WITH DENTURE REMOVAL?

Now you can give them a solution with the CHAIRSIDE® Denture Removal Tool. This easy to use tool is designed to be used by the patient, or a caregiver, to assist with the dislodgement of a dental prosthesis including dentures and implant overdentures.

PATIENT FRIENDLY AND EASY TO USE

- Large ergonomic handle for easy, no slip grasping. Bi-directional design for patients with manual dexterity limitations.
- Angled handle facilitates a comfortable approach to the intraoral environment and the extended neck allows the patient to reach prosthetic flanges while keeping the tool at a patient comfortable level (at or below the chin).
- Rounded tip and no sharp edges to prevent soft tissue trauma or prosthesis damage.

ITEM NUMBER	DESCRIPTION	PRICE
09549	Denture Removal Tool (10 pack)	\$49.99



Insert removal edge under flange of denture.



Use handle of tool to pull down to disengage the overdenture.

GIVE YOUR PATIENTS THE CHAIRSIDE DENTURE REMOVAL TOOL TODAY!
CALL 800.262.2310 OR VISIT WWW.ZESTDENT.COM

ZD ZEST DENTAL SOLUTIONS
 ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS
 2875 Loker Avenue East
 Carlsbad, CA 92010
 USA

TEL: 442.244.4835 OR 800.262.2310
 FAX: 760.743.7975
 EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM

DISTRIBUTED BY: