

TurboTempEZ™

Long-Term, Long-Span Provisional Composite

Making Temporary Restorations EZ

Introducing TurboTemp EZ™, a self-cure provisional crown and bridge composite that offers faster, more reliable, and easier temporary restorations.

FEATURES & BENEFITS

- Reaches a higher flexural strength and hardness in a shorter amount of time than leading brands, allowing for easier trimming and polishing.
- Six essential shades with life-like fluorescence and high polishability offer beautiful results.
- 1:1 cartridge makes for more consistent mixing and curing.
- A predictable set time of less than 5 minutes and reliable handling supports the creation of outstanding temporary restorations.



Work Time	Removal	Set Time
0:45 secs	2:30 mins - 3:00 mins	5:00 mins

Developed as a multifunctional, highly-filled material, TurboTemp EZ is suited for long- and short-term provisional restorations.

INDICATIONS


- Bridges
- Long-span bridges
- Crowns & partial crowns
- Inlays & onlays
- Veneers




Nathaniel Lawson, DMD, PhD
Four-unit bridge fabricated with TurboTemp EZ, A2 shade

4X Higher flexural strength and higher hardness in the critical working stages prevents deformation or shrinkage during handling, trimming and polishing.**

FLEXURAL STRENGTH (MPa)




SHORE D HARDNESS (MPa)




TurboTemp EZ offers provisionals with high compressive strength and gloss, which means durable and beautiful smiles for patients.

COMPRESSIVE STRENGTH, 24 HRS (MPa)



GLOSS



TurboTemp EZ™ Product Ordering Information

Mfg. Part Number	Product Description	MSRP*
90-00085	TurboTemp EZ, A1 Refill (82 gm), 10 mixing tips	\$203.99
90-00086	TurboTemp EZ, A2 Refill (82 gm), 10 mixing tips	\$203.99
90-00087	TurboTemp EZ, A3/3.5 Refill (82 gm), 10 mixing tips	\$203.99
90-00088	TurboTemp EZ, B1 Refill (82 gm), 10 mixing tips	\$203.99
90-00089	TurboTemp EZ, C3 Refill (82 gm), 10 mixing tips	\$203.99
90-00090	TurboTemp EZ, Bleach White, Refill (82 gm), 10 mixing tips	\$203.99
94746-15	Mixing Tips Refill, 15 short, blue mixing tips	\$22.00

*Manufacturer suggested retail price only, pricing may vary according to chosen distributor. TurboTemp EZ is a trademark of Zest Dental Solutions. All rights reserved. All other trademarks and registered trademarks are not trademarks of Zest Dental Solutions.

**Zest Dental Solutions: Data on file.



ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

PHONE: 800.827.7940
FAX: 760.444.4881
EMAIL: DANVILLECS@ZESTDENT.COM
WWW.ZESTDENT.COM