



# A PERFECT MATCH AT YOUR FINGERTIPS

**One composite delivers easy shade matching** and the reliability of long-term durability at a value that better serves your bottom line. ZNano™ is a true nanofilled composite that gives you:

- 5 adaptive shades for easy blending
- Lifelike esthetics with high polishability
- Superior strength & stain resistance
- Ideal nonstick handling for simplified sculpting & placement
- Outstanding effectiveness with fewer steps & less materials
- **Lower price, more savings**



**FREE  
SAMPLE**


REQUEST YOUR FREE SAMPLE TODAY AND DISCOVER  
THE COMPOSITE THAT MATCHES THE WAY YOU PRACTICE.  
**GO TO <https://bit.ly/Try-ZNano>**



## INDICATIONS

ZNano™ Shade-Adaptive Universal Nanocomposite is a highly-filled, light-cured, resin-based dental composite material that is indicated for:

- Direct restorations in cavity Classes I to V
- Composite inlays and onlays
- Reconstruction of damaged teeth, both anterior and posterior
- Diastema closure and correction of the tooth's shape and shade
- Repairing veneers and provisional crown-and-bridge materials



**ZNano™ at 20,000x magnification**

Zirconia-silica nano fillers offer improved shade-matching ability, gloss, and handling.

## Physical Property Chart\*\*

| Brand                                    | Flexural Strength (MPa) | Flexural Modulus (GPa) | Compressive Strength (MPa) | Radiopacity (%Al) | Stain Resistance ( $\Delta E$ , 7 Days Red Wine) | Gloss (GU) | Wear Resistance ( $\mu m$ ) Loss of Thickness |
|--|-------------------------|------------------------|----------------------------|-------------------|--|------------|---|
| ZNano™<br>Zest Dental Solutions          | 134                     | 9                      | 286                        | 205               | 1.5  | 60         | 77  |
| Filtek™ Supreme Ultra 3M                 | 123                     | 10.2                   | 281                        | 223               | 8  | 41         | 206   |
| Estelite® Sigma Quick<br>Tokuyama Dental | 100                     | 8.2                    | 278                        | 161               | 8.5  | 56         | 186   |

## ZNano™ Product Ordering Information

|                           | A1/B1    | A2/B2    | A3/A3.5/B3 | BL/OP    | C2/C3    |
|---------------------------|----------|----------|------------|----------|----------|
| Unit Dose (20x 0.2g tips) | 90-00118 | 90-00119 | 90-00120   | 90-00121 | 90-00122 |
| Syringe (4g)              | 90-00112 | 90-00113 | 90-00114   | 90-00115 | 90-00116 |
| MSRP*                     | \$79.99  | \$79.99  | \$79.99    | \$79.99  | \$79.99  |

\*Manufacturer suggested retail price only, pricing may vary according to chosen distributor. ZNano is a trademark of Zest Dental Solutions. All rights reserved. All other trademarks and registered trademarks are not trademarks of Zest Dental Solutions. \*\*Zest Dental Solutions: Data on file.