

Caries Finder™

Caries Disclosing Dye

Caries Finder helps retain maximum tooth structure and reduce possibility of leaving carious dentin.

FEATURES:

- Prevent Recurrent Caries
- Retain Maximum Tooth Structure
- Available in Red or Green



Caries Finder Green	10ml	80004
Caries Finder Red	10ml	80010
Caries Finder Red Syringable	1gm	90764



HELPS RETAIN MAXIMUM TOOTH STRUCTURE

By precisely staining only carious dentin or enamel Caries Finder dye ensures conservative caries removal. It reduces the possibility of inadvertently leaving carious dentin or of removing excessive dentin. Caries Finder will not stain the following materials:

- Normal healthy dentin
- Inner carious dentin (uninfected) with solid collagen fibers
- Normal healthy enamel

Since Caries Finder only stains the denatured collagen fibers of infected dentin, healthy and remineralizable dentin is easily differentiated.

QUICK AND EASY APPLICATION



ZD ZEST DENTAL SOLUTIONS®
ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM