

Prelude™

Self-Etch, Total Etch & Dual-Cure Adhesive

THE LATEST IN SELF-ETCH TECHNOLOGY

ONE OF THE FASTEST AND EASIEST TO USE


Prelude also offers high bond strengths to uncut enamel and self/dual-cure composites. Prelude Adhesive, used without Primer, also serves as a Total-Etch material. Nano-filled, ethanol-based formulation is highly stable and does not settle out.

ADVANTAGES

- Fast & reliable
- one 10 second application per bottle
- Use of Prelude Link and self/dual-cure composites converts Prelude to dual-cure adhesive
- High bond strengths to uncut enamel!
- Nano-filled, ethanol based formula is highly stable
- Ultra thin film thickness, less than 5 micron
- High 5 minute & 24 hour (34.7MPa) bond strengths



Shear Bond Strength on Dentin at 5min. in self-etch technique



Fastest Application!

