

MicroEtcher™ CD

Intraoral Sandblaster

NO NEED TO COLD STERILIZE YOUR MICROETCHER ANY MORE!




Now we offer a Microetcher™ that plugs right in to your KaVo connector, The Microetcher™ CD.

The introduction of the original Microetcher intraoral sandblaster changed chairside bonding forever. The advent of the Microetcher II allowed for precision microroughening and featured a removable and autoclavable nozzle. We've made major internal improvements to the Microetcher II in order to afford enhanced performance at lower pressure levels. The newly patented nozzle configuration is slimmer in profile allowing greater visibility during intraoral use, can rotate 360 and is interchangeable with the Microetcher II nozzle options. The Microetcher CD is designed to limit abrasive use, in order to minimize overspill in your operatory.

- Connects right to KaVo handpiece connector
- Fully autoclavable body and nozzle
- Compatible with Microetcher II nozzles
- Effective at pressures as low as 45psi
- 30% reduction in abrasive use
- Ergonomically designed nozzle to improve intraoral vision
- Improves bond strength by 400%
- Precise control for intraoral use
- Rapid bonding prep with failures greatly reduced
- Ideal for orthodontic and other appliances

Sand Blasted
Dental Alloy

800x



Bur Roughened
Dental Alloy

800x

