

MicroProphy II™

Air Polisher

Air polishing, the only effective & conservative means of cleaning pits and fissures for sealants. Air polishing rapidly cleans inter proximally, around orthodontic fixtures and into occlusal grooves.


The MicroProphy II is a highly effective, compact and reliable air polishing system. The slim handpiece is finger button controlled for foot pedal free operation.

Highly rated by major reviewers, the MicroProphy II is now improved. Improvements provide increased cleaning efficiency, greater ease of cleaning and use with either a dry or wet spray.

- Wet or Dry Spray
- Finger Button Control
- Compact and Reliable
- Handpiece Rests in Base
- Fully Rotatable, Autoclavable Tips

Unique to MicroProphy, both water and abrasive are contained in the base unit, requiring only connection to air.

The MicroProphy II easily connects to existing "MicroEtcher" air quick connects. Optional kits allow connection to handpiece lines or to install quick disconnects.



ORDERING INFORMATION

MicroProphy II	201684-00
Quick Disconnect Kits	44010

Nozzles

.023 Tip	23623
.026 Tip	23626

Abrasive:

Sodium Bicarbonate B+	17001
Sodium Bicarbonate B	17011

Note: For wet operation use only B+

MICROPROPHY II KIT INCLUDES:

- MicroProphy Unit
- Handpiece with Nozzle (.026)
- Base
- Instructions