

First Quarter™

Fast Setting Impression Material

First Quarter™ F.S. [fast setting (1-1/2 minutes)] impression material is a super hydrophilic VPS impression taking systems that features a unique 1-1/2 minute setting time. As most impressions are that of single units, this material makes a lot of sense to clinicians who desire a faster setting impression taking system for cases that are single units or challenging due to the presence of excessive blood and saliva.

BENEFITS:

- Accurate Impressions in the Presence of Blood & Saliva
- Captures Fine Detail Subgingivally
- Model Can Be Poured for Up to Two Weeks
- Good Asepsis
- Dimensionally Stable
- Ideal for Single Unit C&B Impressions

FEATURES:

- Super Hydrophilic Formulation
- Very High Tear Strength
- Superior Dimensional Stability
- Can Be Sterilized
- Will not Distort in Solution
- Snap Sets in 1-1/2 Minutes



Super Hydrophilic VPS Impression Material Set in 1-1/2 minutes

	2/pk 2-50ml Ctg/6 Tips	4/pk 4-50ml Ctg/12 Tips	Econo Kit 24-50ml Ctg	Super Econo 96-50ml Ctg
Heavy Body	90913	90809	90810	90811
Light Body	89277	89272	89462	90731
Monophase	89386	89387	80008-01	90729
Monophase Choco-Mint	90901	90900	90902	90903