

# Ice Light™

## Light Transmitting Post

ICELIGHT is the next generation of aesthetic posts from Danville!

ICELIGHT features 70% fill of Zir-conia enriched, parallel glass fibers for exceptional strength. Post design is parallel with a tapered end. Its flexibility is close to dentin so stress is distributed, minimizing the risk of cervical fractures.

ICELIGHT is light transmitting for use with light-cured as well as self/dual-cured composites and adhesives.

Posts are also presilanated, making placement rapid. The length is easily trimmed eliminating the stocking of multiple lengths. One reamer, per size, gives the ideal post clearance and allows excess composite to escape. Further, there is no composite post with higher radiopacity!

ICELIGHT kit feature 20 posts in sizes 1mm, 1.2, 1.4 & 1.6mm with reamers.



### NON-LIGHT TRANSMITTING, ICEPOST, ARE ALSO AVAILABLE

**BOTH ARE VALUE PRICED**

VALUE PRICED KITS	REF #
<b>ICEPOST Intro Kit 20/pk</b>	<b>90079-01</b>
Econo Kit 20/pk	90084-01
1.0 mm 10/pk Refill	90082-01
1.2 mm 10/pk Refill	90083-01
1.4 mm 10/pk Refill	90088-01
1.6 mm 10/pk Refill	90089-01

VALUE PRICED KITS	REF #
<b>ICELIGHT Intro Kit 20/pk</b>	<b>90080-01</b>
1.0 mm 10/pk Refill	90378-01
1.2 mm 10/pk Refill	90380-01
1.4 mm 10/pk Refill	90382-01
1.6 mm 10/pk Refill	90384-01