

...its as easy as
1-2-3

ABC Wedge™

Sectional Matrix System

- Restores natural contour
- Ideal contacts
- Minimizes finishing
- Retentive Tines
- Flat Tines
- Contoured stainless steel Matrices
- Occlusal offset
- Stackable

CONTACT RINGS FEATURES:

- **Converging tines** to make the rings far more retentive than parallel tines
- **Inward and outward tines** for placement on either side of the cavity preparation
- **Occlusal offset** to allow stacking and also clearance over obstacles like rubber dam clamps.



ABC Wedge with Contact Matrix™ is the ultimate sectional matrix system, with three improvements over any other matrix system.

The ABC Wedge makes restoring the wide open box or missing cusp easy when using sectional matrix or trimmed Tofflemeyer bands. It is all too frequent that excessive decay results in fractured cusps. Previously, the only alternative was a Tofflemeyer band with the obtrusive retainer, frequently creating difficulty in restoring a tight, biomimetic contact.

ABCWedge Kit (192 Wedges) Ref # 93488

- X-Small (24R, 24L)
- Small (24R, 24L)
- Medium (24R, 24L)
- Large (24R, 24L)
- Mega Rings (2) Gold PN 91304
- Pliers (1) PN 91298

SINGLE REFILLS (192 Wedges)

- | | |
|-------------------------|-------------|
| Large (96 Left/Right) | Ref # 93589 |
| Medium (96 Left/Right) | Ref # 93588 |
| Small (96 Left/Right) | Ref # 93487 |
| X-Small (96 Left/Right) | Ref # 93486 |

SINGLE REFILLS (192 Wedges)

- | | |
|---------------|-------------|
| Large Right | Ref # 93771 |
| Large Left | Ref # 93810 |
| Medium Right | Ref # 93772 |
| Medium Left | Ref # 93811 |
| Small Right | Ref # 93773 |
| Small Left | Ref # 93812 |
| X-Small Right | Ref # 93774 |
| X-Small Left | Ref # 93813 |



ZD ZEST DENTAL SOLUTIONS®
ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM