

The only matrix Ring you'll ever need!




Danville Materials, the originator of the converging tines and sectional matrix, brings you MEGA V, with a unique, inclined hoop design for superior separation, stacking and field of view. Improve your class II, save finishing time and improve your bottom line.



SECTIONAL MATRIX SYSTEM

MEGA V represents the innovation behind the converging tines technology and the superior clinical performance only a sectional matrix system can provide. The use of heat-treated stainless steel, ensures the intended force is delivered to separate the adjacent tooth, freeing the clinician to perform class II restorations, quickly and efficiently. The MEGA V sectional matrix system offers a single ring size intended for molars and premolars, which can be stacked to achieve multiple restorations at one time.

- The inclined hoop design allows for optimal separating force, easy stacking and better field of view
- Durable Nylon® tines provide grip, securely separating while holding the matrix band in place
- MEGA V's slot dimension is maximized to contour to the tooth while supporting the use of multiple wedges
- Matrix bands available in stiff, thin and subgingival to optimize recreation of the perfect margins and contacts
- Compatible with other systems' instruments and wedges



KITS AVAILABLE:

Clinical Kit Plus -Mega V Ring (4 pk), REF 94271	Matrices Thin Flex Small (100 pk) Matrices Thin Flex Large (100 pk) Contact Wedge/ring Plier, Megagrip Matrix Forceps Mega Wedges - Sm, Med, Lg 25@	REF 89432 REF 89434 REF 91298 REF 90486
Clinical Kit - REF 94272	Mega V Ring (4 pk), Matrices Thin Flex Small (100 pk) Matrices Thin Flex Large (100 pk) Mega Wedges - Sm, Med, Lg 25@	REF 89432 REF 89434
Trial Kit - REF 94270	Mega V Ring (2 pk), Matrices Thin Flex Small (50 pk) Matrices Thin Flex Large (50 pk) Mega Wedges - Sm, Med, Lg 25@	REF 90067 REF 90068
Refill Kit - REF 94269	Mega V Ring (2 pk)	

