

Z-Bond™

Primer

Z-Bond is a brush on single bottle primer for zirconia and alumina ceramics used prior to dental filling material for direct restoration (repair) and cements for indirect restoration. While primarily used for zirconia and alumina, the product also provides improved adhesion when used to prepare precious and non-precious metals. Z-Bond is contraindicated on lithium disilicates where silanes perform better.



Z-Bond
5ml
Ref# 93410

ADHESIVE INTRAORAL REPAIR OF DEFECTS

in zirconia and alumina restorations. Conditioning of zirconia and alumina restorations prior to cementation, e.g., crowns, inlays, onlays, etc.