

Danville **MACRO CAB+**

Self-Contained Dust Cabinet

Compact, yet larger, unit offering an efficient and comfortable design where you only insert your fingers with handpiece to work. Bright, LED light provides ideal illumination. Featuring built-in dust collection/filtration system. It is easy to clean and the filter lasts for years. It is quiet, made of heavy gauge aluminum with a baked enamel finish.



MacroCab+
7.9" (h) x 12.5" (w) x 10.5" (d)
Ref # 91556
Disposable Windows
Ref # 93960
Filter
Ref # 91557

- New LED lights with 30,000 hour life that save energy!
- A 50% higher flow fan for better suction power!
- New window design for a better fit - less powder leakage!
- Universal power supply 100-240, 50-60 Hz with UL, CE and TUV approvals.

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance.

94065 Rev A