


# Accolade™

## Thixotropic Flowable Composite

### ADVANTAGES

- Adapts readily, yet does not slump
- 65% fill for superior wear
- Superior polish, 0.7u fill
- Ideal handling, non-sticky
- True Vita Shade match
- Fluoride slow release
- Value priced!

Accolade offers all the highly popular features of StarFlow while resisting slumping. This allows the use of thicker, time saving layers. Accolade provides superb handling, making it easy to work. True Vita shade match, without color shift assures predictable results. New Fluorescent fill gives natural appearance under all lighting conditions. Excellent esthetics make Accolade ideal as a final restorative or as a cavity liner.



ACCOLADE STAYS PUT

WHILE

OTHERS RUN

BEFORE



AFTER



### VITA SHADES AVAILABLE:

A-00op, A-0, A-1, A-2, A-3, A-3.5, A-5, B-1, B-2, C-2, C-4, Dentin & Incisal

### PORCELAIN VENEER SHADES AVAILABLE:

Light, Extra Light, White Opaquer, Translucent, Yellow, & Brown