

TurboTemp™ 2

Veneer Temporaries Made at Turbo Speed

It is especially attractive for anterior veneer temporaries due to its esthetic Vita colors and slight opacity.

Turbo Temp™ 2 has up to twice the flexural strength of competitive products, virtually eliminating failures. This unique, syringeable microfilled composite produces esthetically superior provisional restorations, taking less time than a lab temp. With less than 0.4% shrinkage, Turbo Temp 2 is so accurate that when used in conjunction with a good vinyl polysiloxane such as Start VPS, it makes incomparably fast and accurate temporaries, rarely requiring occlusal adjustment.

ADVANTAGES

- Less than 0.44% shrinkage ensures a precise fit
- Bubbles and margin defects are minimized
- Up to 2 times higher flexural strength than competitive products
- Defects easily repaired by direct bonding
- Fast intraoral setting time



VENEER TEMPORARIES!



STARTER KITS

- 1 - 76gm cartridge
- 10 mix tips
- Dispenser gun
- 76gm
- 10 Mix Tips

| Shade | Ref # |
|-------|-------|
| A-1 | 90343 |
| A-2 | 90344 |
| A-3.5 | 90345 |

| Shade | Ref # |
|-------|-------|
| A-1 | 90156 |
| A-2 | 90157 |
| A-3.5 | 90158 |