

AlphaBite™

Bite Registration Material

Alpha Bite is Danville's new bite registration material that is formulated to the specifications of Blu Mousse®¹. It has the same 'standby' properties, yet is more value priced. It is also easy to carve without breaking. Alpha Bite is made from quality fillers and silicones.

Work time is 45 sec, Set is 75 sec.

Available in 50ml cartridges.

¹Blu Mousse is a trademark of Parkell.

Alpha Bite
(2) 50ml + 6 tips
Ref # 93352



Noginate™

VPS Alginate Substitute

Noginate is a value priced VPS material that may be used as an alginate substitute. Ideal for producing orthodontic models, study models, for creating mouth guards and breaching trays.

Unlike alginate, Noginate is very stable, has a high tear strength and may be disinfected. Noginate is so stable that it can be archived instead of heavy stone models.

Work time: 90 sec. Set: 150 sec.

Value priced, available in 50ml cartridges.

Noginate
(2) 50ml + 6 tips
Ref # 93353



ZD ZEST DENTAL SOLUTIONS™
ZEST | DANVILLE MATERIALS | PERIOSCOPY

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM