

ZEST Anchors Introduces a New Angle for Narrow Diameter Implants

Escondido, CA – (April 7, 2015) – ZEST Anchors, the trusted manufacturer of end-to-end solutions in the resilient implant-retained overdenture market, announced today the release of the SATURNO™ Narrow Diameter Implant System featuring an innovative, Angled O-Ball Implant option and patented, pivoting O-Ring technology.

The new SATURNO Implant is a one-piece narrow diameter implant with a 1.8mm O-Ball connection that incorporates many unique features not found in mini o-ball implants currently on the market today. These features include: a 20° Angled O-Ball Implant to address the challenges associated with implant placement in the maxillary ridge; patented, pivoting cap technology providing maximum resiliency while maintaining a minimized profile; an innovative narrow body implant design incorporating progressive threads for increased primary stability; a self-tapping feature for ease of implant insertion; and a proven RBM surface on the entire length of the implant.

Each SATURNO Implant is available in your choice of 20° Angled or Straight O-Ball configurations; in 2.0, 2.4 or 2.9mm diameters; 10, 12, or 14mm lengths; and cuff heights in 2 or 4mm. The value conscious packaging also includes a Restorative Component Vial with the Attachment Housing and Block-Out Spacer for processing the overdenture.

ZEST is offering a SATURNO Narrow Diameter Implant System product special for clinicians who want to start offering this unique solution to their edentulous patients. Please call 800-262-2310 or visit www.zestanchors.com to get started today!

About ZEST Anchors, LLC

ZEST Anchors is a global leader in the manufacturing and distribution of dental solutions for the treatment of edentulous patients. ZEST's impressive 40 year history has led to product solutions that have enhanced the outcomes for both clinicians and patients alike. ZEST Anchors has received worldwide acclaim for pioneering the pivoting and self-aligning design of the LOCATOR® Attachment System and continues its innovative path by consistently introducing new products, giving clinicians the tools to provide world-class overdenture therapy to their patients. The most recent product innovations include: the LOCATOR Overdenture Implant (LODI) System – the first of its kind narrow diameter implant system with a detachable LOCATOR Abutment; the CHAIRSIDE® Product Portfolio - dental tools and materials that provide an end-to-end solution for implant-retained overdenture modification and processing; and most recently, The SATURNO Narrow Diameter Implant System – a unique one-piece narrow diameter o-ball implant that has a 20° O-Ball Implant option and features patented, pivoting O-Ring technology. ZEST Anchors is located in Escondido, California and has global distribution through OEM implant companies, distributor networks, and a domestic retail sales operation.

Media Contact Information:

Lisa Adams | Tel: 760-743-7744 ext. 155 | Email: lisa.adams@zestanchors.com

###